



Basket for laboratory glassware

of spot-welded wire-mesh 8 x 8 x 1 mm.

Last digit:

-11 Wire, polyethylene white

-13 Stainless steel wire, electrolytically polished

-15 Wire, polyamide white

Note:

This basket is also available with partitions dividing it into 2 or 4 compartments.

Price on request.

Ref.-No.	Ref.-No.	Ref.-No.			
PE, white	Stainless	Polyamide	Lenght mm	Width mm	Height mm
520-01-11	520-01-13	520-01-15	100	100	100
520-03-11	520-03-13	520-03-15	120	120	120
520-05-11	520-05-13	520-05-15	140	140	140
520-07-11	520-07-13	520-07-15	150	150	150
520-08-11	520-08-13	520-08-15	160	160	160
520-11-11	520-11-13	520-11-15	180	180	180
520-12-11	520-12-13	520-12-15	200	200	150
520-13-11	520-13-13	520-13-15	200	200	200
520-14-11	520-14-13	520-14-15	250	150	150
520-15-11	520-15-13	520-15-15	300	200	150
520-17-11	520-17-13	520-17-15	300	200	200
520-18-11	520-18-13	520-18-15	400	200	150
520-19-11	520-19-13	520-19-15	250	250	250
PE, white	Stainless	Polyamide	Length mm	Width mm	Height mm
520-20-11	520-20-13	520-20-15	300	300	300
520-21-11	520-21-13	520-21-15	400	200	200