

Survey:

Application:



Narrow-necked



Material	Soda-lime glassw
Maximum service temperature	approx. 100 °C
Maximum difference in temperature	Δ T approx. 30 K
Microwave use	Yes
Use with pressure and vacuum	No
Cleaning of soda-lime glass	- Maximum cleaning temperature 60 °C - not dishwasher-proof

Wide-mouth



behr Labor-Technik GmbH • Spangerstraße 8 • 40599 Düsseldorf/Germany
Tel.: (+49) (0) 211 – 7 48 47 17 • Fax: (+49) (0) 211 – 7 48 47 48
eMail: info@behr-labor.com • Internet: www.behr-labor.com

Subject to errors and technical changes



Sampling and Sample Storage

Glass Bottles with PTFE Inlay



150104-EN-V3_1_001

Precondition for analytical precision: Sampling and sample storage without contamination

Alterations and contaminations as a result of faulty sampling and/or sample storage inevitably affect subsequent analysis and cause faulty results. These faults cannot be corrected afterwards

Contemporary trace analysis is especially sensitive to faulty sampling and/or sample storage.

Modern measuring instruments record even the smallest contamination and alteration.

behrotest® sampling bottles with PTFE coated inlay in the lid deliver safety for the analyst

The sample gets in touch with nothing but glass and PTFE. Thus there is no danger of contamination.

Bottles from brown glass are available for samples that are sensitive to light.

Narrow-necked and wide-mouth bottles from 30 ml to 1,000 ml cover your requirements


behrotest® sampling bottles with PTFE coated inlay in the lid are available as

- narrow-necked bottles for liquid samples. Additionally to being inert, the inlay in the lid makes the bottles absolutely tight.
- wide-mouth bottles for solid and pasty samples.

The Bottle for Your Samples

The inlay in the lid is coated with PTFE.

So the content of the bottle will have contact with nothing but glass and PTFE.



Narrow-necked Bottles, Brown Glass, Lid With PTFE Coated Inlay

Sampling bottles from partly opaque brown glass. PTFE coated inlay in the lid, thus inert to a large extend and especially tight.

Pack of 10

Type	Volume	Ø Thread mm	Stock no.
NB 50 GT	50 ml	28	B00636012
NB 100 GT	100 ml	28	B00636013
NB 250 GT	250 ml	28	B00636014
NB 500 GT	500 ml	28	B00636015
NB 1000 GT	1000 ml	28	B00636016



NB 50 GT



Narrow-necked Bottles, Clear Glass

Narrow-necked Bottles, Clear Glass, Lid With PTFE Coated Inlay

Sampling bottle from clear glass. PTFE coated inlay in the lid, thus inert to a large extend and especially tight.

Pack of 10

Type	Volume	Ø Thread mm	Stock no.
NK 50 GT	50 ml	28	B00636017
NK 100 GT	100 ml	28	B00636018
NK 250 GT	250 ml	28	B00636019
NK 500 GT	500 ml	28	B00636020
NK 1000 GT	1000 ml	28	B00636021



NK 50 GT



Wide-mouth Bottles, Brown Glass

Wide-mouth Bottles, Brown Glass, Lid With PTFE Coated Inlay

Sampling bottles from partly opaque brown glass. PTFE coated inlay in the lid, thus inert to a large extend and especially tight.

Pack of 10

Type	Volume	Ø Thread mm	Stock no.
RB 30 GT	30 ml	32	B00622064
RB 50 GT	50 ml	32	B00622065
RB 100 GT	100 ml	40	B00622066
RB 250 GT	250 ml	55	B00218524
RB 500 GT	500 ml	55	B00218525
RB 1000 GT	1000 ml	68	B00218526



RB 1000 GT



Wide-mouth Bottles, Clear Glass

Wide-mouth Bottles, Clear Glass, Lid With PTFE Coated Inlay

Sampling bottle from clear glass. PTFE coated inlay in the lid, thus inert to a large extend and especially tight.

Pack of 10

Type	Volume	Ø Thread mm	Stock no.
RK 30 GT	30 ml	32	B00622067
RK 50 GT	50 ml	32	B00622068
RK 100 GT	100 ml	40	B00622069
RK 250 GT	250 ml	55	B00218531
RK 500 GT	500 ml	55	B00218532
RK 1000 GT	1000 ml	68	B00218533



RK 1000 GT