


Condensers acc. to Allihn, with GL threads



Condensers acc. to Allihn, with GL threads

Allihn condensers acc. to DIN 12581 and in addition to DIN, for condensing vapours, with NS socket and NS cone, water connection: GL 14 threads (with two threaded plastic olives), made of DURAN® tubing

Socket & cone (NS)	Jacket length (mm)	Cooling surface (cm ²)	Balls (pcs.)	DIN	Pack (qty.)	Cat. No.
14/23	160	100	6	-	1	5 2232 11
19/26	160	100	6	12581	1	5 2232 12
24/29	160	100	6	-	1	5 2232 13
29/32	160	100	6	-	1	5 2232 14
14/23	250	200	5	-	1	5 2233 11
19/26	250	200	5	-	1	5 2233 12
24/29	250	200	5	-	1	5 2233 13
29/32	250	200	5	12581	1	5 2233 14
24/29	400	300	8	-	1	5 2234 13
29/32	400	300	8	-	1	5 2234 14
29/32	500	400	10	-	1	5 2235 14