


Distillation thermometers



Distillation thermometers

Precision thermometers with cone NS 14/23, top end melted round, frosted glass with black scale,

Range (°C)	Scale (°C)	Built-in length (mm)/filling	Pack (qty.)	Cat. No.
0-150	0,5	55/parafin oil	1	5 0650 01
0-150	0,5	70/parafin oil	1	5 0650 02
0-150	0,5	85/parafin oil	1	5 0650 03
0-250	1	55/isoamyl benzoate	1	5 0650 04
0-250	1	70/isoamyl benzoate	1	5 0650 05
0-250	1	85/isoamyl benzoate	1	5 0650 06