



STOPPERS - "cellulose"

- Manufactured from high grade cellulose which can be sterilized up to 200°C.
- These stoppers are economical, labour and time saving and perfectly tight fitting. They are air permeable and actively breathing.
- Excellent sealing for microbiological samples and tissue cultures in test tubes and flasks.

cellulose stoppers	bottom diameter	top diameter	height	pack quantity
051.24.001	7,5 mm	12,5 mm	30,0 mm	10 pieces
051.24.002	9,0 mm	12,5 mm	30,0 mm	10 pieces
051.24.003	11,0 mm	15,5 mm	32,0 mm	10 pieces
051.24.004	11,5 mm	16,5 mm	40,0 mm	10 pieces
051.24.005	12,5 mm	17,0 mm	40,0 mm	10 pieces
051.24.006	15,5 mm	20,0 mm	40,0 mm	10 pieces
051.24.007	15,5 mm	25,0 mm	45,0 mm	10 pieces
051.24.008	20,0 mm	26,5 mm	50,0 mm	10 pieces
051.24.009	24,0 mm	35,0 mm	55,0 mm	10 pieces



STOPPERS - "turn - over flange" - "rubber"

- Manufactured from rubber or silicone with same technical specifications of other stoppers.
- These stoppers guarantee perfect sealing even with poorly finished rims on glass or plastic bottles.
- The top of the stopper has a septum which can be punctured by a hypodermic needle.
- Very secure fit prevents accidental opening of the flask.

rubber stoppers	silicone stoppers	diameter	pack quantity
051.25.001	051.26.001	7,1 mm	10 pieces
051.25.002	051.26.002	10,2 mm	10 pieces
051.25.003	051.26.003	12,7 mm	10 pieces
051.25.004	051.26.004	14,9 mm	10 pieces
051.25.005	051.26.005	15,9 mm	10 pieces
051.25.006	051.26.006	19,4 mm	10 pieces
051.25.007	051.26.007	23,7 mm	10 pieces
051.25.008	051.26.008	30,7 mm	10 pieces



STOPPERS - "polypropylene"

- Manufactured from polypropylene according to NS joints dimensions for secure closure of ground jointed all kinds of laboratory glassware.

clear stoppers	NS size	height	pack quantity
051.08.010	NS 10/19	23 mm	10 pieces
051.08.012	NS 12/21	26 mm	10 pieces
051.08.014	NS 14/23	30 mm	10 pieces
051.08.019	NS 19/26	31 mm	10 pieces
051.08.024	NS 24/29	37 mm	10 pieces
051.08.029	NS 29/32	42 mm	10 pieces