



Made for your laboratory

Providing High Quality Tissue Processing
Reagents manufactured in Europe



At EpreDia, we pride ourselves on providing high quality tissue processing reagents. To better serve our European customers, we have established an EU-based tissue processing reagent manufacturing site with innovative offerings, extended expiration dates and reduced packaging for you and the environment.

EpreDia is committed to streamlining manufacturing practices to reduce waste and increase on-time delivery.

EpreDia Fixatives

EpreDia is now offering high-quality Neutral Buffered Formalin (NBF) Fixative compliant with IVDR REGULATION (EU) 2017/746 regulation for the histology lab in bottles, canisters, and conveniently packaged pre-filled pods for safety and convenience in your laboratory.

Item Number	Description	Order Size	Certification
5706-LP	Fixative, 10% NBF	12 x 1 L	IVD
5735-LP	Fixative, 10% NBF	4 x 5 L	IVD
5710-LP	Fixative, 10% NBF	10 L	IVD
51201-LP-Y	Prefill, 10% NBF 20/10ml - yellow cap	24 x 10 ml	IVD
51601-LP-Y	Prefill, 10% NBF 60/30ml - yellow cap	35 x 30 ml	IVD
531251-LP-Y	Prefill, 10% NBF 125/60ml - yellow cap	24 x 60 ml	IVD
531801-LP-Y	Prefill, 10% NBF 180/80ml - yellow cap	24 x 80 ml	IVD
535001-LP	Prefill, 10% NBF 500/250ml	250 ml	IVD
531001-LP	Prefill, 10% NBF 1L/500ml	500 ml	IVD

At EpreDia, we know how important a label means to your specimen’s integrity throughout the workflow. Our prefilled pots promote superior sample integrity by optimizing a standardized process from patient biopsy acquisition to pathology. Easy identification labels allow you to record the required patient identifiers for proper labeling to meet your workflow needs.

- Patient’s complete first and last name
- Medical record number
- Date and time of specimen collection
- Source of specimen for each container (i.e. biopsy site, tissue type)



2025 Product Update



Eprelia Alcohols and Xylene

Ensure an efficiently running histology laboratory with a trusted partner. Eprelia now offers high-quality tissue processing and histology reagents such as xylene, alcohol, and isopropanol manufactured in the EU for laboratories, all from one supplier.

Item Number	Description	Order Size	Certification
6215-LP	Denatured Absolute Alcohol	4 x 5 L	IVD
6315-LP	95% Dehydrant - Ethanol	4 x 5 L	IVD
6515-LP	70% Dehydrant - Ethanol	4 x 5 L	IVD
9517-LPR	Dehydrant - Isopropanol	4 x 5 L	RUO
9517	Dehydrant - Isopropanol*	4 x 5 L	IVD
6615-LP	Xylene	New Certification 4 x 5 L	IVD

Eprelia
CE-IVD offering
of ethanols are
manufactured
with a 4 year
shelf life!

Eprelia is minimizing the environmental impact by the reduction of packaging with our products. The **NEW Eprelia 5 liter tissue processing reagents** uses less plastic compared to most standard bottles, reducing waste and sourced in EU.

5 Liter Tissue Processing Reagents

- Unique bottle shape allows more reagent capacity in standard flammable cabinets
- 5-liter fill reduces tissue processor underfill errors that can cause delays in patient diagnosis
- 5-liter taps offered specifically to minimize spills and promote a *safer reagent exchange* for laboratory staff

Eprelia Syntri Safeguard™ Reagents

Syntri Safeguard™ reagents are specifically designed for use on the Revos rotational tissue processor. The innovative built-in RFID technology tracks and logs reagent usage for innovative reagent management, providing a total solution.

Item Number	Description	Order Size	Certification
6215SS-LP	Syntri Safeguard™ Denatured Absolute Alcohol	4 x 5 L	IVD
6615SS	Syntri Safeguard Xylene*	4 x 5 L	IVD
9517SS	Syntri Safeguard Isopropyl Alcohol*	4 x 5 L	IVD
6764515SS	Syntri Safeguard Xylene Substitute*	4 x 5 L	IVD
5735SS	Syntri Safeguard Fixative, 10% NBF*	4 x 5 L	IVD

Using Syntri Safeguard™ RFID technology partnered with Eprelia Revos™ ensures that the correct reagent is exchanged and creates an exchange log to track every bottle loaded, protocol executed, a record of lot number and expiration date – helping to reduce downtime, meet your record keeping needs and facilitate your workload.



*Note: This product is manufactured in the USA.

PINK-Fix for easy visualization of sample

Once tissue samples are removed from the body, starting the fixation process immediately is necessary, specifically when working with small biopsies. Rose-colored formalin in prefilled containers is designed to visualize small, tissue fragments and biopsy-sized specimens better.

- The colored solution improves the visibility of grey or opaque tissue fragments, aiding grossing technicians and pathologists throughout the grossing process.
- Allowing technicians to visualize small biopsies at the bench may provide better efficiency and mitigates the error of tissue loss during processing and embedding.
- Further visualization of small tissue biopsies throughout the workflow process further maintains the integrity of the tissue.



Item Number	Description	Order Size	Certification
5735P-LP	Fixative, 10% PINK-Fix	4 x 5 L	IVD
51201P-LP-Y	Prefill, 10% PINK-Fix 20/10 ml, yellow cap	24 x 10 ml	IVD
51601P-LP-Y	Prefill, 10% PINK-Fix 60/30 ml, yellow cap	35 x 30 ml	IVD
531251P-LP-Y	Prefill, 10% PINK-Fix 125/60 ml, yellow cap	24 x 60 ml	IVD

Another step in standardization of downstream testing

Water is an essential component in many reactions. While it seems convenient to use tap water for laboratory use, this might create issues that are difficult to detect. The type of water used can have permanent adverse effects, like sub-optimal pH, reaction-affecting ions, or similar of your buffers. These conditions could affect your staining morphology, RNA/DNA extraction ability, protein evaluation, and histochemical staining.

DI-water is ideal for laboratories performing downstream testing, like molecular tests or IHC – for standardized and reproducible conditions.



Item Number	Description	Order Size	Certification
92905000-LP	Pure de-mineralised (DI) water	4 x 5 L	General lab supply
92910000-LP	Pure de-mineralised (DI) water	10 L	General lab supply
LP-TAP5L	Dispenser taps adapted to 5L EU cans	EA	General lab supply
LP-TAP10L	Dispenser taps adapted to 10L EU cans	EA	General lab supply

Coming in 2026

Item Number	Description	Order Size	Certification
6115SS-LP	Syntri-Safeguard Xylene	4 x 5 L	IVD
9517-LP	Dehydrant - Isopropanol	4 x 5 L	IVD
9517SS-LP	Syntri Safeguard Dehydrant - Isopropanol	4 x 5 L	IVD

Find out more at [epredia.com](https://www.epredia.com)