

# MICRO-TOUCH™ Sterile

P e r f o r m a n c e   S e r i e s


## Latex sterile medical exam gloves with comfortable inner coating

- **Sterile exam grade comfort:** Sterile latex gloves with a soft, inner polyacrylate coating ensure a comfortable fit and easy donning
- **Secure grip:** Textured fingers provide a safe and secure grip in wet or dry conditions
- **Enhanced sensitivity:** Thin design offers excellent tactile sensitivity and dexterity
- **Reduced skin irritation:** Low residual protein content to reduce risk of contact dermatitis
- **Minimized contamination risk:** Packaged in pairs with a protective double wrapping

\*CAUTION: This product contains natural rubber latex, which may cause allergic reactions

### Recommended For

- Sterile procedures outside of the O.R.
- Standard, moderate risk examination procedures
- Damp donning



### MICRO-TOUCH™ Everyday, Everywhere Protection

MICRO-TOUCH™ premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

# MICRO-TOUCH™ Sterile

Performance Series

## PRODUCT INFORMATION

|                          |                                       |                                         |                    |
|--------------------------|---------------------------------------|-----------------------------------------|--------------------|
| Material                 | Natural Rubber Latex                  | Not Made From Natural Rubber Latex      | No                 |
| Color                    | White                                 | Cuff Length                             | Standard           |
| External Glove Surface   | Fully Textured                        | Freedom from Holes (Inspection level I) | 1.5 AQL            |
| Available Sizes          | S (6.5 - 7), M (7.5 - 8), L (8.5 - 9) | Tested For Use With Chemotherapy Drugs  | No                 |
| Sterile                  | No                                    | Antistatic                              | Not Tested         |
| Glove Length (mm/inches) | 245 / 9.65                            | Product Segmentation                    | Performance Series |

## ORDERING INFORMATION

| Size         | S (6.5 - 7) | M (7.5 - 8) | L (8.5 - 9) |
|--------------|-------------|-------------|-------------|
| Product Code | 553342      | 553343      | 553344      |

## PACKAGING AND STORAGE

|                      |                                                                                                          |
|----------------------|----------------------------------------------------------------------------------------------------------|
| Packaging            | 1 pair per pouch 50 pouches per dispenser 4 dispensers per case 200 gloves per case                      |
| Shelf Life           | 5 Years                                                                                                  |
| Storage Instructions | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition. |

## STANDARDS AND CERTIFICATIONS

CE 2797, Medical Device Regulation (EU) 2017/745 Class I, ISO 13485, ISO 9001, ISO 10993-10:2004, UKCA 0086

## Vulcanization Chemical Accelerators

|                                                                                                                                            |                                    |
|--------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|
| Zinc Diethyldithiocarbamate (ZDEC)                                                                                                         | Zinc Dibutyldithiocarbamate (ZDBC) |
| Only a very small number of users may be sensitive to this ingredient(s) and hence may develop irritant and/or allergic contact reactions. |                                    |

Contact your Ansell representative for ordering or more information.

North America  
US Tel: 800 952 9916  
CA Tel: 1-844-494-7854  
Email:  
insidesalesus@ansell.com

Latin America & Caribbean  
Region  
BR +55 (11) 3356-3100  
latam.medical@ansell.com  
MX: +52 442 296 2050  
latam.medical@ansell.com  
CO: +57 1 288 3247  
latam.medical@ansell.com

Europe, Middle East & Africa  
Tel: +32 (0) 2 528 74 00  
Email: info@ansell.eu

Asia Pacific  
Tel: +603 8310 6688  
Email:  
apac.medical@ansell.com

Australia & New Zealand  
Tel: +61 3 9270 7270  
Email:  
protection@ap.ansell.com