

Wide-necked bottle | easy to fill and clean, compact form, unbreakable



Sturdy bottle with wide-necked opening and screw cap. Particularly suitable for viscous media, powder, granules and pastes. Leakproof.

- ▶ LDPE or PP
- ▶ LDPE coloured brown for light-sensitive liquids

Individual customer-specific printing possible - please ask us!

Wide-necked bottle

Material	Capacity ml	Colour	Dia. mm	Height mm	Neck dia. inside mm	Thread dia. mm	Item no.	
LDPE	50	transparent	40	80	24	32	0318-0050	N3
LDPE	100	transparent	48	95	24	32	0318-0100	N3
LDPE	200	transparent	59	114	30	40	0318-0200	N3
LDPE	250	transparent	64	127	31	40	0318-0250	N3
LDPE	500	transparent	77	158	39	50	0318-0500	N3
LDPE	750	transparent	88	176	39	50	0318-0750	N3
LDPE	1000	transparent	94	207	50	65	0318-1000	N3
LDPE	1500	transparent	108	232	50	65	0318-1500	N3
LDPE	2000	transparent	120	252	50	65	0318-2000	N3
PP	100	transparent	48	95	24	32	0319-0100	N3
PP	250	transparent	64	127	31	40	0319-0250	N3
PP	500	transparent	77	158	39	50	0319-0500	N3
PP	1000	transparent	94	207	50	65	0319-1000	N3
LDPE	250	brown	64	127	31	40	0320-0250	N3
LDPE	500	brown	77	158	39	50	0320-0500	N3
LDPE	1000	brown	94	207	50	65	0320-1000	N3