



Two-way valves




Three-way valves



Two and three-way valves are used for transfer of liquids or gases in order to release the flow, to block it or to change its direction.

- ▶ PVDF: excellent chemical resistance, transparent, suitable for food
- ▶ PP/PE: for industrial use, red/white
- ▶ Operating temperature: PVDF +0 up to +40 °C, PP/PE +5 up to +40 °C
- ▶ Resistant to compression max. 1 bar
- ▶ With mounting plate for easy installation



Two-way valves

Material	For dia. mm	Nominal width mm	Item no.	
PP/PE	5-7	4	8605-0060	N4
PP/PE	7-9	6	8605-0080	N4
PP/PE	9-11	8	8605-0100	N4
PP/PE	11-13	8	8605-0120	N4
PVDF	5-7	4	8607-0060	N4
PVDF	7-9	6	8607-0080	N4
PVDF	9-11	8	8607-0100	N4
PVDF	11-13	8	8607-0120	N4

Three-way valves

Material	For dia. mm	Nominal width mm	Item no.	
PP/PE	5-7	4	8606-0060	N4
PP/PE	7-9	6	8606-0080	N4
PP/PE	9-11	8	8606-0100	N4
PP/PE	11-13	8	8606-0120	N4
PVDF	5-7	4	8608-0060	N4
PVDF	7-9	6	8608-0080	N4
PVDF	9-11	8	8608-0100	N4
PVDF	11-13	8	8608-0120	N4