

Magnetic stirring bars



Cylindrical stirring bars with smooth surface, many uses.

- ▶ PTFE opaque
- ▶ Max. use temperature +270 °C
- ▶ High performance permanent magnets and long lifespan

Magnetic stirring bar set consists of 10 of the most standard bars, complete with crystal clear box.

- ▶ 1 of each 15 x 6 mm, 30 x 7mm
- ▶ 2 of each 20 x 7 mm, 25 x 6 mm, 40 x 7 mm, 50 x 7 mm

Magnetic stirring bars			
Length mm	Dia. mm	Item no.	
15	6	8271-1506 N3	
20	6	8271-2006 N3	
20	7	8271-2007 N3	
25	6	8271-2506 N3	
30	6	8271-3006 N3	
30	7	8271-3007 N3	
40	7	8271-4007 N3	
40	8	8271-4008 N3	
50	7	8271-5007 N3	
50	8	8271-5008 N3	
80	10	8271-8010 N3	

Magnetic stirring bar set			
Model	Item no.		
set incl. ten bars	8271-1000		N3

Stirring paddle



Stirring paddle made of polypropylene for gentle stirring in plastic and glass containers.

- ▶ PP
- ▶ Shatterproof
- ▶ Warp resistant
- ▶ Total length 1,17 m or 1,71 m



Stirring paddle				
Length mm	Width mm	Total length cm	Item no.	
155	130	117	8203-2100	N3
190	167	171	8203-3150	N3