



SIEVE ACCESSORIES

RETSCH OFFERS A WIDE SELECTION OF ACCESSORIES FOR PERFECT SIEVE ANALYSES



SAMPLE DIVIDER

They divide all pourable solids so accurately that the characteristic composition of each fraction of the sample corresponds exactly to that of the original bulk sample.

BENEFITS

- | Dust free and exact dividing
- | Solutions for 1 to 10 samples for batch processing
- | Simple handling



FLUID BED DRYER TG 200

The fluid bed dryer TG 200 permits the gentle drying of organic, inorganic, chemical or pharmaceutical bulk materials.

BENEFITS

- | Short drying times
- | gentle drying, dispersing and mixing also of temperature sensitive materials
- | no localized overheating



ULTRASONIC BATHS UR

The ultrasonic bath range UR includes three sizes for cleaning test sieves and grinding tools quickly and easily.

BENEFITS

- | fast, gentle, and highly efficient cleaning
- | increases working life of test sieves
- | easy to use, saves time and cost



SIEVE ACCESSORIES

RETSCH's large selection of accessories makes the perfect sieve analysis complete.

OPTIONS

- | Test sieve racks
- | Accessories for test sieves, wet sieving and sieving aids
- | Sample dividers ...and much more.

OTHER USEFUL TOOLS

- | Accessories for test sieves
Collecting pans, intermediate pans, intermediate rings, sieve lids and o-rings.
- | Accessories for wet sieving
Collecting pans with outlet, venting rings and o-rings.
- | Sieving aids
Chain rings, agate, rubber or steatite balls, brushes, polyurethane cubes.

www.retsched.com/sieves

ORDER DATA

SIEVING AIDS

32.365.0001		Chain ring for test sieves 200 mm and 203 mm Ø to support horizontal sieving
32.050.0001		Brushes, 3 pieces
32.902.0001		Cubes of polyurethane, 12 x 12 x 12 mm, 10 pieces
32.902.0002		Cubes of polyurethane, 20 x 20 x 20 mm, 10 pieces
32.354.0001		Balls of rubber, 20 mm Ø, 5 pieces
32.354.0002		Balls of agate, 10 mm Ø, 10 pieces
32.354.0004		Balls of steatite, 6 mm Ø, 150 g

TEST SIEVE RACK

32.012.0001		Test Sieve Rack for 10 Test Sieves Ø 200 mm/8", height 50 mm/25 mm
-------------	---	--

ACCESSORIES FOR TEST SIEVES (PANS, RINGS, LIDS)

FOR TEST SIEVES 200 MM Ø, HEIGHT 50 MM

69.720.0050		Collecting pan	stainless steel	200 mm Ø	height 50 mm
69.220.0050		Intermediate pan	stainless steel	200 mm Ø	height 50 mm
69.121.0050		Intermediate ring	stainless steel	200 mm Ø	height 50 mm
69.520.0051		Sieve lid	stainless steel	200 mm Ø	

69.420.0050		Collecting pan with outlet	stainless steel	200 mm Ø	height 50 mm
69.221.0025		Venting ring for wet sieving	stainless steel	200 mm Ø	height 25 mm
05.114.0174		O-ring for test sieves	200 mm Ø		

FOR TEST SIEVES 200 MM Ø, HEIGHT 25 MM

69.720.0025		Collecting pan, stainless steel, 200 mm Ø, height 25 mm			
69.220.0025		Intermediate pan, stainless steel, 200 mm Ø, height 25 mm			
69.121.0025		Intermediate ring, stainless steel, 200 mm Ø, height 25 mm			
69.520.0051		Sieve lid, stainless steel, 200 mm Ø			
69.420.0050		Collecting pan with outlet, stainless steel, 200 mm Ø, height 50 mm			
69.221.0025		Venting ring for wet sieving, stainless steel, 200 mm Ø, height 25 mm			
05.114.0174		O-ring for test sieves, 200 mm Ø			

FOR TEST SIEVES 203 MM Ø / 8" Ø, HEIGHT 2"

69.720.3050		Collecting pan, stainless steel, 8" Ø, height 2"			
69.220.3050		Intermediate pan, stainless steel, 8" Ø, height 2"			
69.121.3050		Intermediate ring, stainless steel, 8" Ø, height 2"			
69.520.3051		Sieve lid, stainless steel, 8" Ø			
69.420.3050		Collecting pan with outlet, stainless steel, 8" Ø, height 2"			

69.221.3025  Venting ring for wet sieving, stainless steel, 8" Ø, height 1"

05.114.0174 O-ring for test sieves, 8" Ø

FOR TEST SIEVES 203 MM Ø / 8" Ø, HEIGHT 1"

69.720.3025  Collecting pan, stainless steel, 8" Ø, height 1"

69.220.3025  Intermediate pan, stainless steel, 8" Ø, height 1"

69.121.3025  Intermediate ring, stainless steel, 8" Ø, height 1"

69.520.3051  Sieve lid, stainless steel, 8" Ø

69.420.3050  Collecting pan with outlet, stainless steel, 8" Ø, height 2"

69.221.3025  Venting ring for wet sieving, stainless steel, 8" Ø, height 1"

05.114.0174 O-ring for test sieves, 8" Ø

FOR TEST SIEVES 100 MM Ø

60.010.000100  Collecting pan, stainless steel, 100 mm Ø, height 40 mm

60.220.000100 Intermediate pan, stainless steel, 100 mm Ø, height 40 mm

60.935.000100  Intermediate ring, stainless steel, 100 mm Ø, height 40 mm

60.107.000100  Sieve lid, stainless steel, 100 mm Ø

60.010.100100  Collecting pan with outlet, stainless steel, 100 mm Ø, height 40 mm

05.114.0045  O-ring for test sieves, 100 mm Ø

FOR TEST SIEVES 305 MM Ø

60.010.000305		Collecting pan, stainless steel, 305 mm Ø, height 40 mm
69.230.0050		Intermediate pan, stainless steel, 305 mm Ø, height 40 mm
60.935.000305		Intermediate ring, stainless steel, 305 mm Ø, height 40 mm
60.107.000305		Sieve lid, stainless steel, 305 mm Ø
69.430.0050		Collecting pan with outlet, stainless steel, 305 mm Ø, height 56 mm
69.321.0050		Venting ring for wet sieving, stainless steel, 305 mm Ø, height 40 mm
05.114.0047		O-ring for test sieves, 305 mm Ø

FOR TEST SIEVES 400 MM Ø

60.010.000400		Collecting pan, stainless steel, 400 mm Ø, height 65 mm
60.220.000400		Intermediate pan, stainless steel, 400 mm Ø, height 65 mm
60.935.000400		Intermediate ring, stainless steel, 400 mm Ø, height 65 mm
60.107.000400		Sieve lid, stainless steel, 400 mm Ø
69.440.0070		Collecting pan with outlet, stainless steel, 400 mm Ø, height 70 mm
69.421.000400		Venting ring for wet sieving, stainless steel, 400 mm Ø, height 65 mm
05.114.0048		O-ring for test sieves, 400 mm Ø

FOR TEST SIEVES 450 MM Ø

60.010.000450		Collecting pan, stainless steel, 450 mm Ø, height 50 mm
---------------	---	---

60.220.000450		Intermediate pan, stainless steel, 450 mm Ø, height 100 mm
60.935.000450		Intermediate ring, stainless steel, 450 mm Ø, height 100 mm
69.545.0050		Sieve lid, stainless steel, 450 mm Ø
69.445.0084		Collecting pan with outlet, stainless steel, 450 mm Ø, height 85 mm
69.421.000450		Venting ring for wet sieving, stainless steel, 450 mm Ø, height 100 mm
05.114.0012		O-ring for test sieves, 450 mm Ø