



Specification

1.16834.0001 TLC Silica gel 60 F₂₅₄ 50 Aluminium sheets 5 x 10 cm

Specifications

Specific surface area (according to BET; 5-Pt. measurement)	480 - 540	m ² /g
Pore volume (N ₂ -isotherm)	0.74 - 0.84	ml/g
d 50 (laser diffraction, size distribution)	9.7 - 11.7	µm
Layer thickness	165 - 235	µm
Deviation of layer thickness per plate	≤ 40	µm

hRf-values

Colour test	
Bleu VIF Organol	11 - 25
Ceres black	34 - 48
Ceres violet BRN	52 - 67

Typical value determined on a conditioned sheet
Eluent: Toluene (45% rel. humidity)

Tom Kupfer

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.