



ISOHEAT®



High Tech Engineering

LabHEAT®

Accessories for heating mantles

The appropriate accessories allow a minimum of effort to ensure a tilt-proof set-up for the heating mantles or to integrate the heating mantles within support walls.

Tripod series KM-DF

for higher stability of the standard heating mantles, the application of the tripod is recommended.


The heating mantle is hereby simply on its four clips in the tripod suspended.



Support ring series KM-TR

to integrate the standard heating mantle within support wall, the application of the support ring is recommended. The heating mantle is able to be suspended on its four clips in the support ring.

With a cross-clamp, the arm of the support ring can be fastened to the vertical or horizontal support wall.



Support clamp type KM-SK

for metal-cased heating mantles and serial heating unit to fix support rods up to Ø13 mm or to integrate the metalcased heating mantles within support wall.



Ordering details

Tripod series KM-DF

Volume ml/l	Ring inside Ø mm	Height mm	Part.-No.
100	110	130	71000100
250	140	150	71000250
500	160	160	71000500
1	188	180	71001000
2	222	210	71002000
3	238	230	71003000
4	261	240	71004000
5	277	240	71005000
6	291	265	71006000
10	334	310	71010000
20	430	390	71020000

Ordering details

Support ring series KM-TR

Volume ml/l	Ring inside Ø mm	Arm lenght mm	Part.-No.
25	90	185	72000025
50	90	185	72000050
100	110	185	72000100
250	140	185	72000250
500	160	185	72000500
1	188	230	72001000
2	222	185	72002000
3	238	230	72003000
4	261	230	72004000

Ordering details

Support clamp KM-SK and Support rods KM-SS

KM-SK	for support rods up to Ø 12 mm	Part.-No. 97000001
KM-SS500	Ø 12 mm x 500 mm, without winding	Part.-No. 97000500
KM-SS1000	Ø 12 mm x 1000 mm, without winding	Part.-No. 97001000



series KM-DF



series KM-TR



type KM-SK with KM-SS



ISOHEAT GmbH
Großer Sand 4
D-76698 Ubstadt-Weiher

Tel.: +49 (0) 72 51 - 93247-0
Fax +49 (0) 72 51 - 93247-29

E-Mail: info@isoheat.de
Internet: www.isoheat.de