



designed for scientists

Montageanleitung / Assembly instructions I-MAG UHC



Lesen Sie die Montageanleitung vor Inbetriebnahme vollständig und beachten Sie die Sicherheitshinweise. Wenn dieses Zubehör mit einem anderen Gerät verwendet wird, beachten Sie auch dessen Bedienungsanleitung.



Read the assembly instructions in its entirety before using the device and follow the safety instructions. If this accessory is used with another device, observe also its operating instructions.



Sicherheitshinweise

Hinweis!

- › Alle Schraubverbindungen müssen sicher angezogen sein.
- › Überprüfen Sie die Montageposition, um ein Umkippen des Gerätes zu vermeiden.

Safety instructions

Notice!

- › All screw connections must be properly tightened.
- › Check the mounting position to avoid the danger of the unit tipping over.

Technische Daten


Werkstoff	Edelstahl 1.4301
Abmessungen (B x T x H)	350 x 266 x 140 mm
Gewicht	3,3 kg

Technische Änderungen vorbehalten!

Technical data


Material	Stainless steel 1.4301
Dimensions (W x D x H)	350 x 266 x 140 mm
Weight	3.3 kg

Subject to technical changes!




Application of I-MAG UHC


① On the desktop





② On the wall





Hinweis: Das Befestigungsmaterial ist den Wandverhältnissen entsprechend auszuwählen und ist nicht im Lieferumfang enthalten.
Note: The fixing screws must be selected according to the wall conditions and will not be included in the scope of delivery.



③ On the vessel




① ② ③ ④



On the rod

④



IKA-Werke GmbH & Co. KG

Janke & Kunkel-Straße 10, 79219 Staufen, Germany • Phone: +49 7633 831-0 • Fax: +49 7633 831-98 • eMail: sales@ika.de

Discover and order the fascinating products of IKA online: www.ika.com



IKAworlwide



IKAworlwide /// #lookattheblue



@IKAworlwide

Technical specifications may be changed without prior notice.
 20000039881_20117357_I-MAG UHC_122022