

[www.saint-tech.lv](http://www.saint-tech.lv)

Do more with OHAUS



SaintTech



# EXPLORER™ Semi-Micro

Built for sophisticated, professional weighing



## Application

Weighing, percent weighing, parts counting, check weighing, dynamic/animal weighing, filling, totalization, formulation, fill weight variation, differential weighing, density determination, peak hold, pipette adjustment, SQC

## Display

Full-colour VGA graphic display (diagonal 5.7"/145mm), detachable, 4-wire resistive touch screen. QWERTY keyboard and numeric keypad to quickly input GLP and GMP data and other application data

## Operation

AC Adapter (included)

## Communication

Easy access communication ports including 2 units of USB port, RS232 and an optional 4th Ethernet port, GLP/GMP data output with real time clock, Direct Data Transfer, print to USB drive

## Construction

Metal base, ABS top housing, stainless steel pan, anti-static glass draftshield with flip-top door, replaceable in-use cover, automatic doors on specific models

## Design Features

AutoCal™ fully-automatic internal adjustment system equipped with 2 internal adjustment weights, routine test and minimum weight determination, fast stabilization time, four touchless sensors, built-in ionizer with automatic doors models, up to 14 operating languages, approved models available, save to USB function, menu lock switch, security bracket, integral weigh below hook for below balance weighing applications, removable stainless steel weighing platform, stability indicator, overload and underload indicators, auto standby



Explorer's weighing cells precision-machined from a solid metal block, together with AutoCal™ provides exceptional accuracy and durability suitable for laboratory and industrial environments.

Features such as four-level user management with password protection, and un-editable system log plus a high level of configurability make Explorer balances well suited for regulated applications.

A detachable terminal with a large colour touchscreen, programmable IR sensors, optional built-in ionizer & automatic doors, and a frameless draftshield make Explorer balances easy & convenient to use.

| Model      | EC Type Approved | Auto Door | Ionizer | Capacity  | Readability (d) | (e) | Platform | Cal. weight, class | Item N°  | Price €  |
|------------|------------------|-----------|---------|-----------|-----------------|-----|----------|--------------------|----------|----------|
| EX125D     |                  |           |         | 52g/120g  | 0,01mg/ 0,1mg   |     | Ø 80mm   | 100g, E2           | 30218998 | 5865,00  |
| EX125      |                  |           |         | 120g      | 0,01mg          |     | Ø 80mm   | 100g, E2           | 30218999 | 6070,00  |
| EX225D     |                  |           |         | 120g/220g | 0,01mg/ 0,1mg   |     | Ø 80mm   | 200g, E2           | 30219000 | 6365,00  |
| EX225D/AD  |                  | •         | •       | 120g/220g | 0,01mg/ 0,1mg   |     | Ø 80mm   | 200g, E2           | 30219001 | 8590,00  |
| EX225/AD   |                  | •         | •       | 220g      | 0,01mg          |     | Ø 80mm   | 200g, E2           | 30219002 | 13710,00 |
| EX125DM    | •                |           |         | 52g/120g  | 0,01mg/ 0,1mg   | 1mg | Ø 80mm   | 100g, E2           | 30139515 | 6305,00  |
| EX125M     | •                |           |         | 120g      | 0,01mg          | 1mg | Ø 80mm   | 100g, E2           | 30139516 | 6520,00  |
| EX225DM    | •                |           |         | 120g/220g | 0,01mg/ 0,1mg   | 1mg | Ø 80mm   | 200g, E2           | 30139517 | 6810,00  |
| EX225DM/AD | •                | •         | •       | 120g/220g | 0,01mg/ 0,1mg   | 1mg | Ø 80mm   | 200g, E2           | 30139518 | 9030,00  |
| EX225M/AD  | •                | •         | •       | 220g      | 0,01mg          | 1mg | Ø 80mm   | 200g, E2           | 30139519 | 14140,00 |

# ADVENTURER™ Semi-Micro

*Intuitive Balances Designed for Routine Weighing*



**EASY OPERATION  
& TIME SAVING**



## Application

Weighing (17 units + custom unit), parts counting, percentage weighing, check weighing, animal/dynamic weighing, totalization/statistics, formulation, density determination, display hold, batch printing

## Display

4.3" (109 mm) Full-colour VGA graphic touchscreen with user-controlled brightness, 6 mechanical buttons

## Operation

AC adapter (included)

## Communication

RS232, USB Host, USB Device (included), GLP/GMP compliant data output with real-time clock

## Construction

Metal base, ABS top housing, stainless steel pan, glass draftshield with two piece top mounted side doors and sliding top door (1mg models only), illuminated up-front level indicator, integral weigh below hook, security bracket, adjustment lock, full housing replaceable in-use cover

## Design Features

Selectable environmental filters, auto tare, user-selectable span adjustment points, software lockout and reset menu, user-selectable communication settings and data print options, user definable project and user IDs, up to 9 operating languages



Adventurer balances feature a color touchscreen, icon-based user interface, and an ergonomic design – making them easy to configure and use.

Features such as specialized weighing modes, multiple connectivity options, and AutoCal™ provide versatility and flexibility for a variety of applications.

Durable construction, large weighing surfaces, a space-saving draftshield design, and full housing in-use cover allow for use in lab, education and industrial environments.

| Model   | EC Type Approved | Capacity | Readability (d) | (e) | Platform | Cal. weight, class | Item N°  | Price € |
|---------|------------------|----------|-----------------|-----|----------|--------------------|----------|---------|
| AX85    |                  | 82g      | 0,01mg          |     | Ø 80mm   | 50g, E2            | 30852308 | 4400,00 |
| AX125D  |                  | 82/120g  | 0,01mg/ 0,1mg   |     | Ø 80mm   | 100g, E2           | 30852307 | 4800,00 |
| AX225D  |                  | 102/220g | 0,01mg/ 0,1mg   |     | Ø 80mm   | 200g, E2           | 30852309 | 5100,00 |
| AX85M   | •                | 85g      | 0,01mg          | 1mg | Ø 80mm   | 50g, E2            | 30852311 | 4660,00 |
| AX125DM | •                | 82/120g  | 0,01mg/ 0,1mg   | 1mg | Ø 80mm   | 100g, E2           | 30852310 | 5090,00 |
| AX225DM | •                | 102/220g | 0,01mg/ 0,1mg   | 1mg | Ø 80mm   | 200g, E2           | 30852312 | 5400,00 |

# PIONEER™ Semi-Micro

Combining Economy and High Performance for Essential Weighing



OPTION

**\$ ECONOMY & ACCURACY**

## Application

Basic weighing, parts counting, percent weighing, animal weighing, density determination, formulation

## Display

Backlit Liquid Crystal Display (LCD); Second line display for additional information or guidance (dot matrix)

## Operation

AC adapter (included)

## Communication

RS232 and USB device (included); GLP/GMP data output with real-time clock

## Construction

Metal base, ABS top housing, removable stainless steel pan, glass draftshield, integrated weigh-below-hook, security bracket, calibration lock and replaceable in-use cover

## Design Features

User-selectable environmental filters and brightness settings, auto-tare, auto-dim, user-selectable span calibration points, software lockout and reset menu, user-selectable communication settings and data print options, user-definable project and user IDs, software overload/underload indicator, stability indicator, four operating languages



The PX offers high accuracy and repeatability for essential weighing applications in laboratory, industrial and education settings at an economical price point.

Featuring a cast metal lower housing, sub-pan and stainless steel weighing pan, the PX is durably constructed for versatile, long-term use.

Pioneer features a second line display for additional information or guidance, a static removal bar for convenient grounding, and USB connectivity.

| Model   | EC Type Approved | AutoCal™ | Capacity | Readability (d) | (e) | Platform | Cal. weight, class | Item N°  | Price € |
|---------|------------------|----------|----------|-----------------|-----|----------|--------------------|----------|---------|
| PX125D  |                  | •        | 52g/120g | 0,01mg/0,1mg    |     | Ø 80mm   | 100g, E2           | 30478196 | 3015,00 |
| PX85    |                  | •        | 82g      | 0,01mg          |     | Ø 80mm   | 50g+20g, E2        | 30485340 | 3105,00 |
| PX225D  |                  | •        | 82g/220g | 0,01mg/0,1mg    |     | Ø 80mm   | 200g, E2           | 30478197 | 3230,00 |
| PX125DM | •                | •        | 52g/120g | 0,01mg/0,1mg    | 1mg | Ø 80mm   | 100g, E2           | 30480168 | 3240,00 |
| PX85M   | •                | •        | 82g      | 0,01mg          | 1mg | Ø 80mm   | 50g+20g, E2        | 30485341 | 3330,00 |
| PX225DM | •                | •        | 82g/220g | 0,01mg/0,1mg    | 1mg | Ø 80mm   | 200g, E2           | 30480169 | 3455,00 |

# Accessories for Semi-Micro Balances

| Accessories                          | Explorer™<br>Semi-Micro | Adventurer™<br>Semi-Micro | Pioneer™<br>Semi-Micro | Item N°  | Price € |
|--------------------------------------|-------------------------|---------------------------|------------------------|----------|---------|
| Density determination kit            | •                       | •                         | •                      | 80253384 | 680,00  |
| Sinker Glass for Density kit         | •                       | •                         | •                      | 83034024 | 355,00  |
| ION-100A EU stand-alone ionizer      | •                       | •                         | •                      | 30095929 | 1460,00 |
| Grid pan                             | •                       |                           |                        | 30298582 | 197,00  |
| Pipette adjustment kit               | •                       |                           |                        | 30252145 | 1230,00 |
| Ethernet interface                   | •                       |                           |                        | 83021082 | 138,00  |
| Interface Kit, RS232-USB             | •                       | •                         | •                      | 30304101 | 78,25   |
| Cable, USB A to Mini USB, Black      | •                       | •                         |                        | 28120263 | 21,20   |
| In-Use Cover for EX                  | •                       |                           |                        | 83033633 | 31,75   |
| In-Use Cover for AX                  |                         | •                         |                        | 30111792 | 33,25   |
| In-Use Cover for PX                  |                         |                           | •                      | 30372546 | 30,00   |
| Dust Cover for EX                    | •                       |                           |                        | 30093334 | 41,00   |
| Anti-theft device (Models over 8 kg) | •                       | •                         | •                      | 80850000 | 68,50   |
| Terminal extension cable 3m          | •                       |                           |                        | 83021083 | 22,60   |
| Terminal extension cable 9m          | •                       |                           |                        | 30078078 | 30,50   |
| SF40A impact printer                 | •                       | •                         | •                      | 30064202 | 715,00  |
| Auxiliary display, AD7-RS            | •                       | •                         | •                      | 30472064 | 96,00   |
| Anti theft device                    | •                       | •                         | •                      | 80850043 | 211,00  |



DENSITY DETERMINATION KIT  
80253384



SINKER GLASS FOR DENSITY KIT  
83034024



ION-100A EU STANDALONE IONIZER  
30095929



GRID PAN  
30298582



PIPETTE ADJUSTMENT KIT  
30252145



ETHERNET INTERFACE  
83021082



INTERFACE KIT, RS232-USB  
30304101



IN-USE COVER  
30372546



IN-USE COVER  
83033633



DUST COVER  
30093334



ANTI THEFT DEVICE  
80850000



TERMINAL EXTENSION CABLE 3M  
83021083



TERMINAL EXTENSION CABLE 9M  
30078078




SF40A IMPACT PRINTER  
30064202



AUXILIARY DISPLAY, AD7-RS  
30472064

# EXPLORER™ Analytical

Highest performance, most technically advanced balance from OHAUS available also with auto door feature



## Application

Weighing (17 units + custom units), parts counting, percentage weighing, animal/dynamic weighing, check weighing, filling, totalization, formulation, differential weighing, density determination, peak hold, ingredient costing, pipette adjustment, fill weight variation, gross/net/tare weighing

## Display

Full-colour VGA graphic display (diagonal 5.7"/145mm), detachable, 4-wire resistive touch screen. QWERTY keyboard and numeric keypad to quickly input GLP and GMP data and other application data

## Operation

AC adapter (included)

## Communication


Easy access communication ports including 2 units of USB port, RS232 and an optional 4th Ethernet port, GLP and GMP including date and time, Direct Data Transfer, print to USB drive

## Construction

Metal base, ABS top housing, stainless steel pan, anti-static glass draftshield with flip-top door, replaceable in-use cover, automatic door on specific models

## Design Features

3-doors frameless draftshield with removable panels with optional automatic draftshield door, modular display design with adjustable viewing angle (attach to base, wall or optional tower), AutoCal™ automatic internal Adjustment, level assist feature with illuminated level indicator, selectable environmental and auto-print settings, integral weigh below hook, protective in-use cover, software lock switch



Equipped with a weighing cell that's precision-machined from a solid metal block, Explorer balances are extremely accurate and durable, and are ideal for laboratory and industrial environments.

Features such as four-level user management with password protection, and un-editable system log plus a high level of configurability make Explorer balances well suited for regulated applications.

A detachable terminal with a large colour touchscreen, programmable IR sensors, automatic doors, and a frameless draftshield make Explorer analytical balances easy and convenient to use.

| Model     | EC Type Approved | Auto Door | Capacity | Readability (d) | (e) | Platform | Cal. weight, class | Item N°  | Price € |
|-----------|------------------|-----------|----------|-----------------|-----|----------|--------------------|----------|---------|
| EX124     |                  |           | 120g     | 0.1mg           |     | Ø 90mm   | 100g, E2           | 83021331 | 3580,00 |
| EX224     |                  |           | 220g     | 0.1mg           |     | Ø 90mm   | 200g, E2           | 83021332 | 3980,00 |
| EX324     |                  |           | 320g     | 0.1mg           |     | Ø 90mm   | 300g, E2           | 83021335 | 4590,00 |
| EX124/AD  |                  | •         | 120g     | 0.1mg           |     | Ø 90mm   | 100g, E2           | 30061976 | 4625,00 |
| EX224/AD  |                  | •         | 220g     | 0.1mg           |     | Ø 90mm   | 200g, E2           | 30061977 | 5015,00 |
| EX324/AD  |                  | •         | 320g     | 0.1mg           |     | Ø 90mm   | 300g, E2           | 30061978 | 5625,00 |
| EX224M    | •                |           | 220g     | 0.1mg           | 1mg | Ø 90mm   | 200g, E2           | 83021333 | 4240,00 |
| EX224M/AD | •                | •         | 220g     | 0.1mg           | 1mg | Ø 90mm   | 200g, E2           | 30061995 | 5280,00 |
| EX324M    | •                |           | 320g     | 0.1mg           | 1mg | Ø 90mm   | 300g, E2           | 83021336 | 4850,00 |
| EX324M/AD | •                | •         | 320g     | 0.1mg           | 1mg | Ø 90mm   | 300g, E2           | 30061996 | 5895,00 |

# ADVENTURER™ Analytical

*Intuitive Balances Designed for Routine Weighing*



**EASY OPERATION  
& TIME SAVING**

## Application

Weighing (17 units + custom unit), parts counting, percentage weighing, check weighing, animal/dynamic weighing, totalization/statistics, formulation, density determination, display hold, batch printing

## Display

4.3" (109 mm) Full-colour VGA graphic touchscreen with user-controlled brightness, 6 mechanical buttons

## Operation

AC adapter (included)

## Communication

RS232, USB Host, USB Device (included), GLP/GMP compliant data output with real-time clock

## Construction

Metal base, ABS top housing, stainless steel pan, glass draftshield with two piece top mounted side doors and sliding top door (1mg models only), illuminated up-front level indicator, integral weigh below hook, security bracket, adjustment lock, full housing replaceable in-use cover

## Design Features

Selectable environmental filters, auto tare, user-selectable span adjustment points, software lockout and reset menu, user-selectable communication settings and data print options, user definable project and user IDs, up to 9 operating languages



Adventurer balances feature a color touchscreen, icon-based user interface, and an ergonomic design – making them easy to configure and use.

Features such as specialized weighing modes, multiple connectivity options, and AutoCal™ provide versatility and flexibility for a variety of applications.

Durable construction, large weighing surfaces, a space-saving draftshield design, and full housing in-use cover allow for use in lab, education and industrial environments.

| Model   | EC Type<br>Approved | AutoCal™ | Capacity | Readability<br>(d) | (e) | Platform | Cal. weight, class | Item N°  | Price € |
|---------|---------------------|----------|----------|--------------------|-----|----------|--------------------|----------|---------|
| AX124   |                     | •        | 120g     | 0.1mg              |     | Ø 90mm   | 100g, E2           | 30122610 | 2240,00 |
| AX224   |                     | •        | 220g     | 0.1mg              |     | Ø 90mm   | 200g, E2           | 30122611 | 2340,00 |
| AX324   |                     | •        | 320g     | 0.1mg              |     | Ø 90mm   | 200g+100g, E2      | 30122613 | 2410,00 |
| AX124M  | •                   | •        | 120g     | 0.1mg              | 1mg | Ø 90mm   | 100g, E2           | 30100620 | 2450,00 |
| AX224M  | •                   | •        | 220g     | 0.1mg              | 1mg | Ø 90mm   | 200g, E2           | 30100621 | 2560,00 |
| AX324M  | •                   | •        | 320g     | 0.1mg              | 1mg | Ø 90mm   | 200g+100g, E2      | 30100622 | 2630,00 |
| AX124/E |                     |          | 120g     | 0.1mg              |     | Ø 90mm   | 100g, E2           | 30122657 | 1970,00 |
| AX224/E |                     |          | 220g     | 0.1mg              |     | Ø 90mm   | 200g, E2           | 30122612 | 2080,00 |

# PIONEER™ Analytical

Combining Economy and High Performance for Essential Weighing



**\$ECONOMY  
& ACCURACY**

## Application

Basic weighing, parts counting, percent weighing, animal weighing, density determination, formulation

## Display

Backlit Liquid Crystal Display (LCD); second line display for additional information or guidance (dot matrix)

## Operation

AC adapter (included)

## Communication

RS232 and USB device (included); GLP/GMP data output with real-time clock

## Construction

Metal base, ABS top housing, removable stainless steel pan, glass draftshield, integrated weigh-below-hook, security bracket, calibration lock and replaceable in-use cover

## Design Features

User-selectable environmental filters and brightness settings, auto-tare, auto-dim, user-selectable span calibration points, software lockout and reset menu, userselectable communication settings and data print options, user-definable project and user IDs, software overload/underload indicator, stability indicator, four operating languages



The PX offers high accuracy and repeatability for essential weighing applications in laboratory, industrial and education settings at an economical price point.

Featuring a cast metal lower housing, sub-pan and stainless steel weighing pan, the PX is durably constructed for versatile, long-term use.

Pioneer features a second line display for additional information or guidance, a static removal bar for convenient grounding, and USB connectivity.

| Model   | EC Type<br>Approved | Internal<br>Adjustment | Capacity | Readability<br>(d) | (e) | Platform | Cal. weight,<br>class | Item N°  | Price € |
|---------|---------------------|------------------------|----------|--------------------|-----|----------|-----------------------|----------|---------|
| PX124   |                     | •                      | 120g     | 0,1mg              |     | Ø 90mm   | 100g, E2              | 30429803 | 1530,00 |
| PX224   |                     | •                      | 220g     | 0,1mg              |     | Ø 90mm   | 200g, E2              | 30429804 | 1640,00 |
| PX124M  | •                   | •                      | 120g     | 0,1mg              | 1mg | Ø 90mm   | 100g, E2              | 30430010 | 1745,00 |
| PX224M  | •                   | •                      | 220g     | 0,1mg              | 1mg | Ø 90mm   | 200g, E2              | 30430013 | 1855,00 |
| PX124/E |                     |                        | 120g     | 0,1mg              |     | Ø 90mm   | 100g, E2              | 30429810 | 1335,00 |
| PX224/E |                     |                        | 220g     | 0,1mg              |     | Ø 90mm   | 200g, E2              | 30429811 | 1450,00 |

# PR SERIES Analytical

*Designed for Routine Weighing Applications in Your Workplace*



**ROUTINE WEIGHING  
FOR ANY BUDGET**

## Application

Basic weighing, parts counting, percent weighing

## Display

Backlit Liquid Crystal Display (LCD)

## Operation

AC adapter (included)

## Communication

RS232 (one port included)

## Construction

ABS top housing, removable stainless steel pan, glass draftshield with sliding top door, integrated weigh-below-hook, security bracket, calibration lock, replaceable in-use cover

## Design features

User-selectable environmental filters and brightness settings, auto-tare, auto-dim, user-selectable span calibration points, software lockout and reset menu, user-selectable communication settings and data print options, user-definable project and user IDs, software overload/underload indicator, stability indicator



Basic Functionality for Routine Weighing Applications.

Designed for Uncomplicated Operation with Easy-to-Use Display and Interface.

Smart Design and Durable Construction.

| Model   | EC Type<br>Approved | Internal<br>Adjustment | Capacity | Readability<br>(d) | (e) | Platform | Cal. weight,<br>class | Item N°  | Price € |
|---------|---------------------|------------------------|----------|--------------------|-----|----------|-----------------------|----------|---------|
| PR124   |                     | •                      | 120g     | 0,1mg              |     | Ø 90mm   | 100g, E2              | 30430067 | 1170,00 |
| PR224   |                     | •                      | 220g     | 0,1mg              |     | Ø 90mm   | 200g, E2              | 30430069 | 1250,00 |
| PR124M  | •                   | •                      | 120g     | 0,1mg              | 1mg | Ø 90mm   | 100g, E2              | 30524384 | 1330,00 |
| PR224M  | •                   | •                      | 220g     | 0,1mg              | 1mg | Ø 90mm   | 200g, E2              | 30524385 | 1430,00 |
| PR124/E |                     |                        | 120g     | 0,1mg              |     | Ø 90mm   | 100g, E2              | 30430068 | 998,00  |
| PR224/E |                     |                        | 220g     | 0,1mg              |     | Ø 90mm   | 200g, E2              | 30430070 | 1075,00 |

# Accessories for Analytical Balances

| Accessories                                              | Explorer™ Analytical | Adventurer™ Analytical | Pioneer™ Analytical | PR Series Analytical | Item N°  | Price € |
|----------------------------------------------------------|----------------------|------------------------|---------------------|----------------------|----------|---------|
| Density determination kit                                | •                    | •                      | •                   |                      | 80253384 | 680,00  |
| Sinker Glass for Density kit                             | •                    | •                      | •                   |                      | 83034024 | 355,00  |
| ION-100A EU stand-alone ionizer                          | •                    | •                      | •                   | •                    | 30095929 | 1460,00 |
| In-Use Cover for EX                                      | •                    |                        |                     |                      | 83033633 | 31,75   |
| Full housing in use cover for 0.1mg and 1mg model for AX |                      | •                      |                     |                      | 30111792 | 33,25   |
| In-Use Cover for PX                                      |                      |                        | •                   |                      | 30372546 | 30,00   |
| In-Use Cover for PR                                      |                      |                        |                     | •                    | 30372547 | 30,00   |
| Dust Cover for EX and AX (0.1mg)                         | •                    | •                      |                     |                      | 30093334 | 41,00   |
| Ethernet interface                                       | •                    |                        |                     |                      | 83021082 | 138,00  |
| Cable, USB A to Mini USB, Black                          | •                    |                        |                     |                      | 28120263 | 21,20   |
| Interface Kit, RS232-USB                                 | •                    | •                      | •                   | •                    | 30304101 | 78,25   |
| Cable, RS422, Extension, 3m                              | •                    |                        |                     |                      | 83021083 | 22,60   |
| Cable, RS422, Extension, 9m                              | •                    |                        |                     |                      | 30078078 | 30,50   |
| Auxiliary display, AD7-RS                                |                      | •                      | •                   | •                    | 30472064 | 96,00   |
| Interface Cable, USB (Type A to B)                       |                      | •                      |                     |                      | 83021085 | 22,60   |
| Anti theft device                                        | •                    | •                      | •                   | •                    | 80850043 | 211,00  |
| SF40A impact printer                                     | •                    | •                      | •                   | •                    | 30064202 | 715,00  |
| Carrying Case, EX                                        | •                    |                        |                     |                      | 30251106 | 630,00  |
| Carrying Case, AX                                        |                      | •                      |                     |                      | 30251107 | 555,00  |
| Carrying Case, PX                                        |                      |                        | •                   |                      | 30477091 | 424,00  |



DENSITY DETERMINATION  
80253384



SINKER GLASS FOR DENSITY KIT  
83034024



ION-100A EU STANDALONE IONIZER  
30095929



IN-USE COVER EX  
83033633



DUST COVER  
30093334



ETHERNET INTERFACE  
83021082



ANTI-THEFT DEVICE  
80850000



INTERFACE KIT, RS232-USB  
30304101



CABLE, RS422, EXTENSION, 3M  
83021083



CABLE, RS422, EXTENSION, 9M  
30078078



AUXILIARY DISPLAY, AD7-RS  
30472064



FULL HOUSING IN USE COVER AX  
30111792



INTERFACE CABLE, USB (TYPE A TO B)  
83021085



IN-USE COVER PX  
30372546



IN-USE COVER PR  
30372547



SF40A IMPACT PRINTER  
ER 30064202

# EXPLORER™ Precision

*Intelligent Performance No Matter How Difficult the Task!*



**INTELLIGENT PERFORMANCE**

## Application

Weighing (17 units + custom units), parts counting, percentage weighing, animal/dynamic weighing, check weighing, filling, totalization, formulation, differential weighing, density determination, peak hold, ingredient costing, pipette adjustment (except for Explorer High Capacity), gross/net/tare weighing

## Display

High resolution 5.7" (145mm diagonal) full colour VGA touch display including QWERTY, detachable and numeric key pads plus weighing capacity indication

## Operation

AC adapter (included) or rechargeable battery (High Capacity models only)

## Communication

Easy access communication ports including 2 units of USB port, RS232 and an optional 4th Ethernet port, GLP and GMP including date and time, Direct Data Transfer, print to USB drive

## Construction

Modular design, metal base, ABS top housing, stainless steel pan, glass draftshield with top mounted side doors and flip/sliding top door (1mg models only), illuminated up-front level indicator, integral weigh below hook, security bracket, Adjustment lock, four touchless sensors, full housing in-use cover

## Design Features

Selectable environmental filters, auto tare, user-selectable span Adjustment points, software lockout and reset menu, user-selectable communication settings and data print options, user-definable project and user IDs, auto standby, up to 14 operating language



Equipped with a weighing cell that's precision-machined from a solid metal block, Explorer balances are extremely accurate and durable, and are ideal for laboratory and industrial environments.

Features such as four-level user management with password protection, and un-editable system log plus a high level of configurability make Explorer balances well suited for regulated applications.

A detachable terminal with a large colour touchscreen, programmable IR sensors, and a frameless draftshield make Explorer precision balances easy and convenient to use.

| Model    | EC Type Approved | AutoCal™ | Capacity | Readability (d) | (e)   | Platform      | Cal. weight, class | Item N°  | Price € |
|----------|------------------|----------|----------|-----------------|-------|---------------|--------------------|----------|---------|
| EX223    |                  | •        | 220g     | 1mg             |       | Ø 130mm       | 200g, F1           | 83021341 | 2590,00 |
| EX423    |                  | •        | 420g     | 1mg             |       | Ø 130mm       | 400g, F1           | 83021342 | 3345,00 |
| EX623    |                  | •        | 620g     | 1mg             |       | Ø 130mm       | 500g, F1           | 83021348 | 3895,00 |
| EX1103   |                  | •        | 1100g    | 1mg             |       | Ø 130mm       | 1000g, F1          | 83021350 | 4865,00 |
| EX2202   |                  | •        | 2200g    | 0.01g           |       | 190mm x 200mm | 2000g, F1          | 83021356 | 2565,00 |
| EX4202   |                  | •        | 4200g    | 0.01g           |       | 190mm x 200mm | 4000g, F1          | 83021357 | 3310,00 |
| EX6202   |                  | •        | 6200g    | 0.01g           |       | 190mm x 200mm | 5000g, F1          | 83021360 | 3885,00 |
| EX10202  |                  | •        | 10200g   | 0.01g           |       | 190mm x 200mm | 10000g, F1         | 83021362 | 4545,00 |
| EX6201   |                  | •        | 6200g    | 0.1g            |       | 190mm x 200mm | 5000g, F1          | 83021361 | 2440,00 |
| EX10201  |                  | •        | 10200g   | 0.1g            |       | 190mm x 200mm | 10000g, F1         | 83021365 | 3025,00 |
| EX423M   | •                | •        | 420g     | 1mg             | 0.01g | Ø 130mm       | 400g, F1           | 83021343 | 3575,00 |
| EX1103M  | •                | •        | 1100g    | 1mg             | 0.01g | Ø 130mm       | 1000g, F1          | 83021351 | 5095,00 |
| EX4202M  | •                | •        | 4200g    | 0.01g           | 0.1g  | 190mm x 200mm | 4000g, F1          | 83021358 | 3535,00 |
| EX10202M | •                | •        | 10200g   | 0.01g           | 0.1g  | 190mm x 200mm | 10000g, F1         | 83021363 | 4770,00 |
| EX10201M | •                | •        | 10200g   | 0.1g            | 0.1g  | 190mm x 200mm | 10000g, F1         | 83021366 | 3250,00 |
| EX223/E  |                  |          | 220g     | 1mg             |       | Ø 130mm       | 200g, F1           | 30061979 | 2290,00 |
| EX423/E  |                  |          | 420g     | 1mg             |       | Ø 130mm       | 400g, F1           | 30061980 | 3045,00 |
| EX2202/E |                  |          | 2200g    | 0.01g           |       | 190mm x 200mm | 2000g, F1          | 30061981 | 2265,00 |
| EX4202/E |                  |          | 4200g    | 0.01g           |       | 190mm x 200mm | 4000g, F1          | 30061982 | 3015,00 |
| EX6202/E |                  |          | 6200g    | 0.01g           |       | 190mm x 200mm | 5000g, F1          | 30061983 | 3585,00 |
| EX6201/E |                  |          | 6200g    | 0.1g            |       | 190mm x 200mm | 5000g, F1          | 30061984 | 2150,00 |
| EX12001  | •                | •        | 12000g   | 0.1g            |       | 377mm x 311mm | 10kg, F1           | 30078072 | 4345,00 |
| EX24001  | •                | •        | 24000g   | 0.1g            |       | 377mm x 311mm | 20kg, F1           | 30078073 | 4935,00 |
| EX35001  | •                | •        | 35000g   | 0.1g            |       | 377mm x 311mm | 10+20kg, F1        | 30078074 | 5980,00 |
| EX12001M | •                | •        | 12000g   | 0.1g            | 1g    | 377mm x 311mm | 10kg, F1           | 30057067 | 4575,00 |
| EX24001M | •                | •        | 24000g   | 0.1g            | 1g    | 377mm x 311mm | 20kg, F1           | 30057068 | 5165,00 |
| EX35001M | •                | •        | 35000g   | 0.1g            | 1g    | 377mm x 311mm | 10+20kg, F1        | 30057069 | 6220,00 |

High Capacity

# ADVENTURER™ Precision

## Intuitive Balances Designed for Routine Weighing



**EASY OPERATION  
& TIME SAVING**

### Application

Weighing (17 units + custom unit), parts counting, percentage weighing, check weighing, animal/dynamic weighing, totalization/statistics, formulation, density determination, display hold, batch printing

### Display

4.3" (109 mm) Full-colour VGA graphic touchscreen with user-controlled brightness, 6 mechanical buttons

### Operation

AC adapter (included)

### Communication

RS232, USB Host, USB Device (included), GLP/GMP compliant data output with real-time clock

### Construction

Metal base, ABS top housing, stainless steel pan, glass draftshield with two piece top mounted side doors and sliding top door (1mg models only), illuminated up-front level indicator, integral weigh below hook, security bracket, adjustment lock, full housing replaceable in-use cover

### Design Features

Selectable environmental filters, auto tare, user-selectable span adjustment points, software lockout and reset menu, user-selectable communication settings and data print options, user definable project and user IDs, up to 9 operating languages



Adventurer balances feature a color touchscreen, icon-based user interface, and an ergonomic design – making them easy to configure and use.

Features such as specialized weighing modes, multiple connectivity options, and AutoCal™ provide versatility and flexibility for a variety of applications.

Durable construction, large weighing surfaces, a space-saving draftshield design, and full housing in-use cover allow for use in lab, education and industrial environments.

| Model     | EC Type Approved | AutoCal™ | Capacity | Readability (d) | (e)   | Platform      | Cal. weight, class | Item N°  | Price € |
|-----------|------------------|----------|----------|-----------------|-------|---------------|--------------------|----------|---------|
| AX223     | •                | •        | 220g     | 1mg             |       | Ø 130mm       | 200g, F1           | 30127960 | 1670,00 |
| AX423     | •                | •        | 420g     | 1mg             |       | Ø 130mm       | 200g+200g, F1      | 30127962 | 1830,00 |
| AX623     | •                | •        | 620g     | 1mg             |       | Ø 130mm       | 500g, F1           | 30805890 | 1970,00 |
| AX422     | •                | •        | 420g     | 0,01g           |       | 175mm x 195mm | 200g+200g, F1      | 30206079 | 1350,00 |
| AX822     | •                | •        | 820g     | 0,01g           |       | 175mm x 195mm | 500g+100g, F1      | 30206081 | 1510,00 |
| AX1502    | •                | •        | 1520g    | 0,01g           |       | 175mm x 195mm | 1000g+500g, F1     | 30127966 | 1550,00 |
| AX2202    | •                | •        | 2200g    | 0,01g           |       | 175mm x 195mm | 2000g, F1          | 30127968 | 1670,00 |
| AX4202    | •                | •        | 4200g    | 0,01g           |       | 175mm x 195mm | 2000g+2000g, F1    | 30122616 | 1830,00 |
| AX6202    | •                | •        | 6200g    | 0,01g           |       | 175mm x 195mm | 5000g, F1          | 30805892 | 1970,00 |
| AX2201    | •                | •        | 2200g    | 0,1g            |       | 175mm x 195mm | 2000g, F1          | 30206083 | 1300,00 |
| AX4201    | •                | •        | 4200g    | 0,1g            |       | 175mm x 195mm | 2000g+2000g, F1    | 30127970 | 1530,00 |
| AX8201    | •                | •        | 8200g    | 0,1g            |       | 175mm x 195mm | 5000g+1000g, F1    | 30122619 | 1630,00 |
| AX223M    | •                | •        | 220g     | 1mg             | 0,01g | Ø 130mm       | 200g, F1           | 30100623 | 1880,00 |
| AX423M    | •                | •        | 420g     | 1mg             | 0,01g | Ø 130mm       | 200g+200g, F1      | 30100624 | 2050,00 |
| AX623M    | •                | •        | 620g     | 1mg             | 0,01g | Ø 130mm       | 500g, F1           | 30805896 | 2190,00 |
| AX1502M   | •                | •        | 1520g    | 0,01g           | 0,1g  | 175mm x 195mm | 1000g+500g, F1     | 30100626 | 1770,00 |
| AX2202M   | •                | •        | 2200g    | 0,01g           | 0,1g  | 175mm x 195mm | 2000g, F1          | 30100627 | 1880,00 |
| AX4202M   | •                | •        | 4200g    | 0,01g           | 0,1g  | 175mm x 195mm | 2000g+2000g, F1    | 30100628 | 2050,00 |
| AX6202M   | •                | •        | 6200g    | 0,01g           | 0,1g  | 175mm x 195mm | 5000g, F1          | 30805897 | 2200,00 |
| AX8201M   | •                | •        | 8200g    | 0,1g            | 1g    | 175mm x 195mm | 5000g+1000g, F1    | 30122651 | 1860,00 |
| AX223/E   |                  |          | 220g     | 1mg             |       | Ø 130mm       | 200g, F1           | 30127961 | 1360,00 |
| AX423/E   |                  |          | 420g     | 1mg             |       | Ø 130mm       | 200g+200g, F1      | 30127963 | 1530,00 |
| AX623/E   |                  |          | 620g     | 1mg             |       | Ø 130mm       | 500g, F1           | 30805891 | 1670,00 |
| AX422/E   |                  |          | 420g     | 0,01g           |       | 175mm x 195mm | 200g+200g, F1      | 30206080 | 1040,00 |
| AX822/E   |                  |          | 820g     | 0,01g           |       | 175mm x 195mm | 500g+100g, F1      | 30206082 | 1210,00 |
| AX1502/E  |                  |          | 1520g    | 0,01g           |       | 175mm x 195mm | 1000g+500g, F1     | 30127967 | 1230,00 |
| AX2202/E  |                  |          | 2200g    | 0,01g           |       | 175mm x 195mm | 2000g, F1          | 30127969 | 1360,00 |
| AX4202/E  |                  |          | 4200g    | 0,01g           |       | 175mm x 195mm | 2000g+2000g, F1    | 30122617 | 1530,00 |
| AX6202/E  |                  |          | 6200g    | 0,01g           |       | 175mm x 195mm | 5000g, F1          | 30805893 | 1670,00 |
| AX2201/E  |                  |          | 2200g    | 0,1g            |       | 175mm x 195mm | 2000g, F1          | 30206084 | 977,00  |
| AX4201/E  |                  |          | 4200g    | 0,1g            |       | 175mm x 195mm | 2000g+2000g, F1    | 30127971 | 1210,00 |
| AX8201/E  |                  |          | 8200g    | 0,1g            |       | 175mm x 195mm | 5000g+1000g, F1    | 30122620 | 1340,00 |
| AX12001   |                  | •        | 12000g   | 0,1g            |       | 175mm x 195mm | 10000g, F1         | 30805894 | 1710,00 |
| AX12001M  | •                | •        | 12000g   | 0,1g            | 1g    | 175mm x 195mm | 10000g, F1         | 30805898 | 1950,00 |
| AX12001/E |                  |          | 12000g   | 0,1g            |       | 175mm x 195mm | 10000g, F1         | 30805895 | 1400,00 |

High Capacity

# PIONEER™ Precision

Combining Economy and High Performance for Essential Weighing



**ECONOMY & ACCURACY**

## Application

Basic weighing, parts counting, percent weighing, animal weighing, density determination, formulation

## Display

Backlit Liquid Crystal Display (LCD); Second line display for additional information or guidance (dot matrix)

## Operation

AC adapter (included)

## Communication

RS232 and USB device (included); GLP/GMP data output with real-time clock

## Construction

Metal base, ABS top housing, removable stainless steel pan, glass draftshield (1 mg models only), integrated weigh-below-hook, security bracket, calibration lock and replaceable in-use cover

## Design Features

User-selectable environmental filters and brightness settings, auto-tare, auto-dim, user-selectable span calibration points, software lockout and reset menu, user-selectable communication settings and data print options, user-definable project and user IDs, software overload/underload indicator, stability indicator, four operating languages



The PX offers high accuracy and repeatability for essential weighing applications in laboratory, industrial and education settings at an economical price point.

Featuring a cast metal lower housing, sub-pan and stainless steel weighing pan, the PX is durably constructed for versatile, long-term use.

Pioneer features a second line display for additional information or guidance, a static removal bar for convenient grounding, and USB connectivity.

| Model      | EC Type Approved | Internal Adjustment | Capacity | Readability (d) | (e)   | Platform      | Cal. weight, class | Item N°  | Price € |
|------------|------------------|---------------------|----------|-----------------|-------|---------------|--------------------|----------|---------|
| PX323      |                  | •                   | 320g     | 1mg             |       | Ø 120mm       | 200g+100g, F1      | 30429805 | 1145,00 |
| PX623      |                  | •                   | 620g     | 1mg             |       | Ø 120mm       | 500g+100g, F1      | 30619383 | 1280,00 |
| PX3202     |                  | •                   | 3200g    | 0,01g           |       | Ø 180mm       | 2kg+1kg, F1        | 30429807 | 1145,00 |
| PX6202     |                  | •                   | 6200g    | 0,01g           |       | Ø 180mm       | 5kg+1kg, F1        | 30619385 | 1260,00 |
| PX4201     |                  | •                   | 4200g    | 0,1g            |       | Ø 180mm       | 2kg+2kg, F1        | 30429809 | 1110,00 |
| PX323M     | •                | •                   | 320g     | 1mg             | 0,01g | Ø 120mm       | 200g+100g, F1      | 30430015 | 1315,00 |
| PX623M     | •                | •                   | 620g     | 1mg             | 0,01g | Ø 120mm       | 500g+100g, F1      | 30619386 | 1460,00 |
| PX3202M    | •                | •                   | 3200g    | 0,01g           | 0,1g  | Ø 180mm       | 2kg+1kg, F1        | 30430021 | 1350,00 |
| PX6202M    | •                | •                   | 6200g    | 0,01g           | 0,1g  | Ø 180mm       | 5kg+1kg, F1        | 30714491 | 1465,00 |
| PX4201M    | •                | •                   | 4200g    | 0,1g            | 0,1g  | Ø 180mm       | 2kg+2kg, F1        | 30430024 | 1280,00 |
| PX323/E    |                  |                     | 320g     | 1mg             |       | Ø 120mm       | 200g+100g, F1      | 30429812 | 920,00  |
| PX623/E    |                  |                     | 620g     | 1mg             |       | Ø 120mm       | 500g+100g, F1      | 30619382 | 1050,00 |
| PX3202/E   |                  |                     | 3200g    | 0,01g           |       | Ø 180mm       | 2kg+1kg, F1        | 30429814 | 930,00  |
| PX6202/E   |                  |                     | 6200g    | 0,01g           |       | Ø 180mm       | 5kg+1kg, F1        | 30619384 | 1175,00 |
| PX4201/E   |                  |                     | 4200g    | 0,1g            |       | Ø 180mm       | 2kg+2kg, F1        | 30429816 | 875,00  |
| PX120001   |                  | •                   | 12000g   | 0,1g            |       | 193mm x 248mm | 10kg, F1           | 30868072 | 1220,00 |
| PX120001/E |                  |                     | 12000g   | 0,1g            |       | 193mm x 248mm | 10kg, F1           | 30868073 | 960,00  |

High Capacity

# PR SERIES Precision

*Designed for Routine Weighing Applications in Your Workplace*



**ROUTINE WEIGHING  
FOR ANY BUDGET**



## Application

Basic weighing, parts counting, percent weighing

## Display

Backlit Liquid Crystal Display (LCD)

## Operation

AC adapter (included)

## Communication

RS232 (one port included)

## Construction

ABS top housing, removable stainless steel pan, glass draftshield (1mg models only) with sliding top door, integrated weigh-below-hook, security bracket, calibration lock

## Design features

User-selectable environmental filters and brightness settings, auto-tare, auto-dim, user-selectable span calibration points, software lockout and reset menu, user-selectable communication settings and data print options, user-definable project and user IDs, software overload/underload indicator, stability indicator

Basic Functionality for Routine Weighing Applications.

Designed for Uncomplicated Operation with Easy-to-Use Display and Interface.

Smart Design and Durable Construction.

| Model    | EC Type Approved | Internal Adjustment | Capacity | Readability (d) | (e)   | Platform | Cal. weight, class | Item N°  | Price € |
|----------|------------------|---------------------|----------|-----------------|-------|----------|--------------------|----------|---------|
| PR223    |                  | •                   | 220g     | 1mg             |       | Ø 120mm  | 200g, F1           | 30492730 | 869,00  |
| PR423    |                  | •                   | 420g     | 1mg             |       | Ø 120mm  | 200g+200g, F1      | 30492731 | 879,00  |
| PR523    |                  | •                   | 520g     | 1mg             |       | Ø 120mm  | 500g, F1           | 30492732 | 931,00  |
| PR1602   |                  | •                   | 1600g    | 0,01g           |       | Ø 180mm  | 1kg+500g, F1       | 30492733 | 885,00  |
| PR2202   |                  | •                   | 2200g    | 0,01g           |       | Ø 180mm  | 2kg, F1            | 30492734 | 895,00  |
| PR4202   |                  | •                   | 4200g    | 0,01g           |       | Ø 180mm  | 2kg+2kg, F1        | 30492735 | 946,00  |
| PR4201   |                  | •                   | 4200g    | 0,1g            |       | Ø 180mm  | 2kg+2kg, F1        | 30492736 | 833,00  |
| PR6201   |                  | •                   | 6200g    | 0,1g            |       | Ø 180mm  | 5kg+1kg, F1        | 30498176 | 854,00  |
| PR223M   | •                | •                   | 220g     | 1mg             | 0,01g | Ø 120mm  | 200g, F1           | 30524386 | 1020,00 |
| PR423M   | •                | •                   | 420g     | 1mg             | 0,01g | Ø 120mm  | 200g+200g, F1      | 30524387 | 1030,00 |
| PR523M   | •                | •                   | 520g     | 1mg             | 0,01g | Ø 120mm  | 500g, F1           | 30524608 | 1075,00 |
| PR1602M  | •                | •                   | 1600g    | 0,01g           | 0,1g  | Ø 180mm  | 1kg+500g, F1       | 30524609 | 1035,00 |
| PR2202M  | •                | •                   | 2200g    | 0,01g           | 0,1g  | Ø 180mm  | 2kg, F1            | 30524610 | 1055,00 |
| PR4202M  | •                | •                   | 4200g    | 0,01g           | 0,1g  | Ø 180mm  | 2kg+2kg, F1        | 30524611 | 1090,00 |
| PR4201M  | •                | •                   | 4200g    | 0,1g            | 0,1g  | Ø 180mm  | 2kg+2kg, F1        | 30524612 | 972,00  |
| PR6201M  | •                | •                   | 6200g    | 0,1g            | 1g    | Ø 180mm  | 5kg+1kg, F1        | 30524613 | 998,00  |
| PR223/E  |                  |                     | 220g     | 1mg             |       | Ø 120mm  | 200g, F1           | 30430071 | 710,00  |
| PR423/E  |                  |                     | 420g     | 1mg             |       | Ø 120mm  | 200g+200g, F1      | 30430072 | 720,00  |
| PR523/E  |                  |                     | 520g     | 1mg             |       | Ø 120mm  | 500g, F1           | 30430073 | 766,00  |
| PR1602/E |                  |                     | 1600g    | 0,01g           |       | Ø 180mm  | 1kg+500g, F1       | 30430074 | 725,00  |
| PR2202/E |                  |                     | 2200g    | 0,01g           |       | Ø 180mm  | 2kg, F1            | 30430075 | 736,00  |
| PR4202/E |                  |                     | 4200g    | 0,01g           |       | Ø 180mm  | 2kg+2kg, F1        | 30430076 | 782,00  |
| PR4201/E |                  |                     | 4200g    | 0,1g            |       | Ø 180mm  | 2kg+2kg, F1        | 30430077 | 664,00  |
| PR6201/E |                  |                     | 6200g    | 0,1g            |       | Ø 180mm  | 5kg+1kg, F1        | 30430078 | 684,00  |

# Accessories for Precision Balances

| Accessories                                                   | Explorer™<br>Precision | Explorer™ Precision<br>High Capacity | Adventurer™<br>Precision | Pioneer™<br>Precision | PR Series<br>Precision | Item N°  | Price € |
|---------------------------------------------------------------|------------------------|--------------------------------------|--------------------------|-----------------------|------------------------|----------|---------|
| Density determination kit (for 1mg only)                      | •                      |                                      | •                        | •                     |                        | 80253384 | 680,00  |
| Sinker Glass for Density kit                                  | •                      |                                      | •                        | •                     |                        | 83034024 | 355,00  |
| Auxiliary display, AD7-RS                                     | •                      | •                                    | •                        | •                     | •                      | 30472064 | 96,00   |
| In-Use Cover for EX                                           | •                      | •                                    |                          |                       |                        | 83033633 | 31,75   |
| Full housing in-use cover for AX (1mg models)                 |                        |                                      | •                        |                       |                        | 30111792 | 33,25   |
| Full housing in-use cover for AX (0.01g and 0.1g models)      |                        |                                      | •                        |                       |                        | 30111777 | 33,25   |
| In-Use Cover for PX (0.01g and 0.1g models)                   |                        |                                      |                          | •                     |                        | 30759721 | 30,00   |
| In-Use Cover for PX                                           |                        |                                      |                          | •                     |                        | 30372546 | 30,00   |
| In-Use Cover for PR                                           |                        |                                      |                          |                       | •                      | 30372547 | 30,00   |
| Ethernet interface                                            | •                      | •                                    |                          |                       |                        | 83021082 | 138,00  |
| Draft shield kit for EX 0,01g and 0,1g models                 | •                      |                                      |                          |                       |                        | 83021084 | 334,00  |
| Tower Kit                                                     | •                      |                                      |                          |                       |                        | 83021102 | 86,75   |
| Tower Kit                                                     |                        | •                                    |                          |                       |                        | 30078082 | 86,75   |
| Anti theft device (models over 12kg)                          |                        | •                                    |                          |                       |                        | 80850000 | 68,50   |
| Anti theft device                                             | •                      |                                      | •                        | •                     | •                      | 80850043 | 211,00  |
| Cable, USB A to Mini USB, Black                               | •                      |                                      |                          |                       |                        | 28120263 | 21,20   |
| Interface Kit, RS232-USB                                      | •                      | •                                    | •                        | •                     | •                      | 30304101 | 78,25   |
| Cable, RS422, Extension, 3m                                   | •                      | •                                    |                          |                       |                        | 83021083 | 22,60   |
| Cable, RS422, Extension, 9m                                   | •                      | •                                    |                          |                       |                        | 30078078 | 30,50   |
| SF40A impact printer                                          | •                      | •                                    | •                        | •                     | •                      | 30064202 | 715,00  |
| Battery Kit                                                   |                        | •                                    |                          |                       |                        | 30041295 | 222,00  |
| Adjustable rolling feet                                       |                        | •                                    |                          |                       |                        | 30041470 | 320,00  |
| Interface Cable, USB (Type A to B)                            |                        |                                      | •                        |                       |                        | 83021085 | 22,60   |
| Carrying Case, AX                                             |                        |                                      | •                        |                       |                        | 30251108 | 485,00  |
| Carrying Case, PX                                             |                        |                                      |                          | •                     |                        | 30477092 | 384,00  |
| Pan, 12kg, 300 x 225mm, PX<br>(only for High Capacity models) |                        |                                      |                          | •                     |                        | 30868077 | 35,00   |



DENSITY DETERMINATION KIT  
80253384



SINKER GLASS FOR DENSITY KIT  
83034024



AUXILIARY DISPLAY, AD7-RS  
30472064



IN-USE COVER  
83033633



ETHERNET INTERFACE  
83021082



DRAFT SHIELD KIT  
83021084



TOWER KIT  
83021102



ANTI THEFT DEVICE  
80850000



INTERFACE CABLE, USB (TYPE A TO B)  
83021085



INTERFACE KIT, RS232-USB  
30304101



CABLE, RS422, EXTENSION, 3M  
83021083



FULL HOUSING IN USE COVER  
30111777



SF40A IMPACT PRINTER  
ER30064202



IN-USE COVER  
30372546



BATTERY KIT  
30041295



ADJUSTABLE ROLLING FEET  
30041470

## Placing an order couldn't be easier

### Online Ordering via our OHAUS Partner Shop

Our user-friendly Partner Shop makes ordering and tracking as simple and convenient as possible.

- ✓ Orders are received directly into our warehouse and production site for speed and efficiency
- ✓ Once your order is placed you will receive instant order acknowledgment and a confirmed delivery date within 24 hours so that you can provide your customers with accurate information

### Ordering via fax, e-mail or phone to our Sales Support Team

- ✓ Once your order is placed you will receive a confirmed delivery date within 24 hours via e-mail.

## Track your order easily

### Online order tracking

- ✓ If you place your order via our Partner Shop you can also track it simply and efficiently online and you can also view your previous orders

### Order tracking via our Sales Support Team

- ✓ Our Sales Support Team can tell you the status of your order via phone, fax or e-mail

## Terms & Conditions

Our general terms and conditions of sale can be found within our Partner Shop. Please visit [www.ohaus.com](http://www.ohaus.com) and select your country. Standard delivery is free for orders over € 400.

Below € 400 a handling fee of € 22 will be charged.

## Contact

### Technical Support Centre

e-mail: [tsc@ohaus.com](mailto:tsc@ohaus.com)

Tel: 0041 22 567 53 20

Fax: 0041 22 567 53 21

### Sales Support Centre

e-mail: [ssc@ohaus.com](mailto:ssc@ohaus.com)

Tel: 0041 22 567 53 19

Fax: 0041 22 567 53 21

OHAUS Europe GmbH  
Heuwinkelstrasse 3,  
8606 Nänikon,  
Switzerland

[www.ohaus.com](http://www.ohaus.com)



Saint Tech

[www.saint-tech.lv](http://www.saint-tech.lv)



[www.ohaus.com](http://www.ohaus.com)