

**KERN
316-07**

Date: 20.02.2025

ENGLISH

Individual weight, OIML Class E2 100 g, Compact form, stainless steel polished

CONTACT Address

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

 Phone and fax

Phone : +49 - [0] 7433 - 9933 - 0
Fax : +49 - [0] 7433 - 9933 - 149

 Internet

Email : info@kern-sohn.com
Web : www.kern-sohn.com

DESCRIPTION

- Individual weight, OIML Class E2 100 g, Compact form, stainless steel polished (OIML)
- For this weight, we strongly recommend a box due to the very accurate weight class You will find a selection in the accessories section

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	316-7
Model-Serie	316
Assortment	KERN
Product Type	Test weights
Product group	OIML E2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761094895
Customs tariff number	84239090

Technical data - Primary

Nominal value	100 g
Tolerance (OIML)	0.16 mg
Case aluminium	317-70-600
Case plastic	317-70-400
Case wood	317-70-100

Technical data - Construction

Material	stainless steel polished
----------	--------------------------

Technical data - Measuring system

Nominal value	100 g
---------------	-------

OIML class E2

OIML conform yes

Tolerance (OIML) 0.16 mg

Density (kg/m³) 7950

Volume (cm³) 12.5

Technical data - Packing & Shipping

Dimensions packaging (W×D×H) 22×22×35 mm

Gross weight 0.10 kg

Net weight 0.10 kg

Shipping method Parcel service

Packaging component - by weight - plastics 0 g

Delivery time 1 d

Technical data - Verification

Verification possible no

Optional Services

DAkkS Calibration 962-337

FUNCTIONS

Standard



ACCESSORIES

Model	Description
317-070-100	Wooden Box for Set of Weights KERN 317-070-100
317-070-400	Plastic Weight Box KERN 317-070-400
317-070-600	Aluminium Protective Box KERN 317-070-600

SERVICES

Model	Description
921-337	Guaranteed Plus-deviation (? +1/3 OIML) KERN 921-337
962-337R	DAkkS Recalibration KERN 962-337R
920-33	Marking by etching KERN 920-33
962-337	DAkkS Calibration KERN 962-337
920-331	Marking by etching (last digits of the serial number) KERN 920-331

RELATED ARTICLES DATA

Model name	Nominal value	Tolerance (OIML)	Case aluminium	Case plastic	Case wood
316-01	1 g	0.03 mg	317-10-600	317-20-400	317-10-100
316-02	2 g	0.04 mg	317-20-600	317-20-400	317-20-100
316-03	5 g	0.05 mg	317-30-600	317-30-400	317-30-100
316-04	10 g	0.06 mg	317-40-600	317-40-400	317-40-100
316-05	20 g	0.08 mg	317-50-600	317-50-400	317-50-100
316-06	50 g	0.1 mg	317-60-600	317-60-400	317-60-100
316-14	10 kg	16 mg	317-140-600	317-140-400	317-140-100
316-07	100 g	0.16 mg	317-70-600	317-70-400	317-70-100
316-08	200 g	0.3 mg	317-80-600	317-80-400	317-80-100
316-09	500 g	0.8 mg	317-90-600	317-90-400	317-90-100
316-11	1 kg	1.6 mg	317-110-600	317-110-400	317-110-100
316-12	2 kg	3 mg	317-120-600	317-120-400	317-120-100
316-13	5 kg	8 mg	317-130-600	317-130-400	317-130-100