

PRODUCT BROCHURE

KERN 337-03

Date: 28.12.2025

ENGLISH



Individual weight, OIML Class F2, 5 g, Knob, stainless steel fine turned

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- Individual weight, OIML Class F2, 5 g, Knob, stainless steelfine turned (OIML)
- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	337-03
Model-Serie	337
Assortment	KERN
Product Type	Test weights
Product group	OIML F2 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761006324
Customs tariff number	84239090

Technical data - Primary

Nominal value	5 g
Tolerance (OIML)	0,5 mg
Case aluminium	317-030-600
Case plastic	347-030-400
Case wood	337-030-200

Technical data - Construction

Material	stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	5 g
OIML class	F2
OIML conform	yes
Tolerance (OIML)	0,5 mg
Density (kg/m ³)	7950
Volume (cm ³)	0,625

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	8×8×16 mm
Gross weight	0.05 kg
Net weight	0.05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	no
-----------------------	----

Optional Services

DAkkS Calibration	962-433
-------------------	---------

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>337-030-200</u>	Wooden Box for Set of Weights KERN 337-030-200
<u>347-030-400</u>	Plastic Weight Box KERN 347-030-400
<u>317-030-600</u>	Aluminium Protective Box KERN 317-030-600
<u>315-242</u>	Tweezers KERN 315-242

SERVICES

Model	Description
<u>337-030-200</u>	Wooden Box for Set of Weights KERN 337-030-200
<u>347-030-400</u>	Plastic Weight Box KERN 347-030-400
<u>317-030-600</u>	Aluminium Protective Box KERN 317-030-600
<u>315-242</u>	Tweezers KERN 315-242

RELATED ARTICLES DATA

Model name	Nominal value	Tolerance (OIML)	Case aluminium	Case plastic	Case wood
<u>337-01</u>	1 g	0,3 mg	317-010-600	347-030-400	337-010-200
<u>337-02</u>	2 g	0,4 mg	317-020-600	347-030-400	337-020-200
<u>337-03</u>	5 g	0,5 mg	317-030-600	347-030-400	337-030-200
<u>337-04</u>	10 g	0,6 mg	317-040-600	347-050-400	337-040-200
<u>337-05</u>	20 g	0,8 mg	317-050-600	347-050-400	337-050-200
<u>337-06</u>	50 g	1 mg	317-060-600	347-070-400	337-060-200

Model name	Nominal value	Tolerance (OIML)	Case aluminium	Case plastic	Case wood
<u>337-07</u>	100 g	1,6 mg	317-070-600	347-070-400	337-070-200
<u>337-08</u>	200 g	3 mg	317-080-600	347-080-400	337-080-200
<u>337-09</u>	500 g	8 mg	317-090-600	347-090-400	337-090-200
<u>337-11</u>	1 kg	16 mg	317-110-600	347-110-400	337-110-200
<u>337-12</u>	2 kg	30 mg	317-120-600	347-120-400	337-120-200
<u>337-13</u>	5 kg	80 mg	317-130-600	347-130-400	337-130-200
<u>337-14</u>	10 kg	160 mg	317-140-600	347-140-400	337-140-200
<u>337-141</u>	10 kg	160 mg	not available	not available	337-141-200
<u>337-15</u>	20 kg	300 mg	317-150-600	-	337-150-200
<u>337-151</u>	20 kg	300 mg	not available	not available	337-151-200
<u>337-16</u>	50 kg	800 mg	317-160-600	-	337-160-200
<u>337-161</u>	50 kg	800 mg	not available	not available	337-161-200
<u>337-131</u>	5 kg	80 mg	-	-	337-131-200