

PRODUCT BROCHURE

KERN 347-04

Date: 28.12.2025

ENGLISH



Individual weight, OIML Class M1, 10 g, Knob, stainless steel fine turned

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

DESCRIPTION

- Individual weight, OIML Class M1, 10 g, Knob, stainless steelfine turned (OIML)
- For the verification you need a box. Please add a box for single weights.

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	347-04
Model-Serie	347
Assortment	KERN
Product Type	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	404576 1008670
Customs tariff number	84239090

Technical data - Primary

Nominal value	10 g
Tolerance (OIML)	2 mg
Case aluminium	317-040-600
Case plastic	347-050-400
Case wood	337-040-200

Technical data - Construction

Material	stainless steel fine turned
----------	-----------------------------

Technical data - Measuring system

Nominal value	10 g
OIML class	M1
OIML conform	yes
Tolerance (OIML)	2 mg
Density (kg/m ³)	7950
Volume (cm ³)	1,25

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	10×10×19 mm
Gross weight	0.05 kg
Net weight	0.05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification

Verification possible	yes
-----------------------	-----

Optional Services

DAkkS Calibration	962-634
-------------------	---------

FUNCTIONS

Standard



Option



ACCESSORIES

Model	Description
<u>337-040-200</u>	Wooden Box for Set of Weights KERN 337-040-200
<u>347-050-400</u>	Plastic Weight Box KERN 347-050-400
<u>317-040-600</u>	Aluminium Protective Box KERN 317-040-600
<u>315-242</u>	Tweezers KERN 315-242

SERVICES

Model	Description
<u>337-040-200</u>	Wooden Box for Set of Weights KERN 337-040-200
<u>347-050-400</u>	Plastic Weight Box KERN 347-050-400
<u>317-040-600</u>	Aluminium Protective Box KERN 317-040-600
<u>315-242</u>	Tweezers KERN 315-242

RELATED ARTICLES DATA

Model name	Nominal value	Tolerance (OIML)	Case aluminium	Case plastic	Case wood
<u>347-01</u>	1 g	1 mg	317-010-600	347-030-400	337-010-200
<u>347-02</u>	2 g	1,2 mg	317-020-600	347-030-400	337-020-200
<u>347-03</u>	5 g	1,6 mg	317-030-600	347-030-400	337-030-200
<u>347-04</u>	10 g	2 mg	317-040-600	347-050-400	337-040-200
<u>347-05</u>	20 g	2,5 mg	317-050-600	347-050-400	337-050-200
<u>347-06</u>	50 g	3 mg	317-060-600	347-070-400	337-060-200
<u>347-07</u>	100 g	5 mg	317-070-600	347-070-400	337-070-200
<u>347-08</u>	200 g	10 mg	317-080-600	347-080-400	337-080-200
<u>347-09</u>	500 g	25 mg	317-090-600	347-090-400	337-090-200
<u>347-11</u>	1 kg	50 mg	317-110-600	347-110-400	337-110-200

Model name	Nominal value	Tolerance (OIML)	Case aluminium	Case plastic	Case wood
<u>347-12</u>	2 kg	100 mg	317-120-600	347-120-400	337-120-200
<u>347-13</u>	5 kg	250 mg	317-130-600	347-130-400	337-130-200
<u>347-14</u>	10 kg	500 mg	317-140-600	347-140-400	337-140-200
<u>347-445-100</u>	10 g	2 mg	-	-	-
<u>347-141</u>	10 kg	500 mg	not available	-	337-141-200
<u>347-151</u>	20 kg	1 g	not available	-	337-151-200
<u>347-161</u>	50 kg	2,5 g	not available	-	337-161-200
<u>347-445-100-N00</u>	-	2 mg	-	-	-
<u>347-022-005</u>	20 kg	1 g	-	-	-
<u>347-042-008</u>	40 kg	1,5 g	-	-	-
<u>347-052-010</u>	50 kg	2,5 g	-	-	-
<u>347-062-022</u>	60 kg	3 g	-	-	-
<u>347-131</u>	5 kg	250 mg	-	-	337-131-200