

Unicare Unitrile (GS005*)

blue · medical grade · heavy weight · nitrile/vinyl hybrid examination gloves

Physical Specifications:

Material: Nitrile and vinyl
 (latex free)
 Colour: Blue
 Length: 240 mm
 Weight: 5.5 g
 Thickness: 0.08 mm
 Finish: Powder free
 Fit: Ambidextrous
 Sterile: No

Order Codes and Size Information:

GS0051 | Extra small
 GS0052 | Small
 GS0053 | Medium
 GS0054 | Large
 GS0055 | Extra large

Packaging:

Each dispenser box contains 100 gloves
 Each case contains 10 dispenser boxes

Classification:



AQL 1.5 (Pinholes)



EN455 (Medical Grade)



Food Safe
 (not for use with fatty foods)

Features

- AQL 1.5
- Free from latex
- Free from powder
- Tested to EN1186
- Manufactured using a blend of high quality vinyl and nitrile
- Beaded cuff
- Heavy-weight construction

Benefits

- Meets medical standards for pinholes
- Suitable for those with latex allergies (Type I sensitisation)
- Reduces risks of latex allergies (Type I sensitisation) and minimises particulate contamination
- Suitable for contact with non-fatty food types
- The unique blend of vinyl and nitrile makes the glove more comfortable and provides better dexterity than traditional vinyl gloves
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Heavy-weight construction is tough, durable and long lasting

Sectors

- Care home work
- Cleaning and janitorial
- Dentistry
- Dental hygiene
- Emergency services
- Food processing
- General assembly
- Healthcare
- Hospital and clinical duties
- Janitorial
- Local Authority work
- Painting and decorating
- Printing



For more about the Unicare range, please click on the image



UNIGLOVES UK
 Registered in England No: 4010200
 VAT registration 791817787

3 Ambley Green, Gillingham Business Park,
 Gillingham, ME8 0NJ, UK
 +44 (0)800 049 6602
 enquiries@unigloves.co.uk
 unigloves.co.uk

UNIGLOVES GERMANY
 Arzt- & Klinikbedarf
 Handelsgesellschaft mbH

Camp-Spich-Straße 71
 D-53842 Troisdorf-Spich,
 +49 (0) 2241-9323-0
 unigloves@unigloves.de
 unigloves.de

