

KERN OZL466T241



The flexible, affordable all-rounder with zoom function for schools, training companies, inspection authorities and laboratories



Construction

Eye distance (max):	75
Eye distance (min):	55
Optical system:	Greenough
Rotation of tube 360°:	yes
Tube:	Trinocular
Tube angle:	45°

Functions

Supported operating system:	Android
-----------------------------	---------

Power supply

Input voltage:	100 - 240 V
Mains adapter/adaptor included:	EURO

Category

Brand:	Optics
--------	--------

Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	590 x 440 x 390 mm
Gross weight:	7,050 kg
Net weight:	5,075 kg
Shipping method:	Parcel service