


Sterile Polychloroprene Cleanroom Gloves

BioClean Ultimate Sterile Polychloroprene Gloves offer good tactility for precision work, and comfort for prolonged use. Offering resistance to a range of chemicals the Ultimate is also latex-free and powder-free.

- Powder-free & latex-free
- Tested against ASTM D6978 for cytotoxic drugs
- ESD properties
- Beaded cuff for strength
- Chemical resistant
- Double-donnable
- Non-particulating EasyTear packaging



Industries

- Controlled and Critical Environments
- Production and Manufacturing
- Lab and R&D
- Sterile Processing Department



BUPS

Sterile Polychloroprene Cleanroom
Gloves

TECHNICAL DATA SHEET

PRODUCT INFORMATION

	BUPS
Material	Neoprene (Polychloroprene)
Color	Natural
Shape	Hand Specific
Cuff	Beaded
Manufacturing/QMS Audit Standards	Manufacturing QMS Audit Standards ISO 9001 PPE Regulation 2016 425 Module D
Regulatory/Standards Compliance	ASTM D6978, CE 0598, EN 420:2003 + A1:2009, EN 421:2010, EN ISO 374-1:2016, EN ISO 374-5:2016, EN 420:2003 + A1:2009, EN 455 Part 2, Food Contact, ISO 9001:2015, Category III
Packaging	One pair per inner PE wallet; one wallet per sealed EasyTear PE pouch; 10 pouches per sealed outer PE bag; 20 outer bags per lined carton (200 pairs)
Storage	Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.
Country of Origin	Malaysia
Available sizes	6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0
Powder Content	Powder-Free
External Glove Surface	Textured Fingers and Palm
Sterilization Method	GAMMA irradiation (25 kGy)
Sterilization Minimum Dose	25kGy
Sterility Assurance Level	10 ⁻⁶
Cleanroom Compatibility	ISO Class 4 & EU GMP Grade A
Shelf Life	Five (5) years from date of manufacture.
Tested for use with Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
Protein Level	N/A: contains no natural rubber latex



BUPS

Sterile Polychloroprene Cleanroom Gloves

PHYSICAL PROPERTIES

	Typical Values							Testing Method
	6	65	7	75	8	85	9	EN420
Length (mm/in)	300/12							
Palm Width	77/3	83/3.2	89/3.5	95/3.7	102/4	108/4.2	115/4.5	
Freedom from Holes	0.65 AQLPerformance Level 3							EN 374-2
Typical Particle Count > 0.5µm (counts / cm²)	<1300							IEST-RP-CC005.4
Minimum Single Wall Palm Thickness (mm/mil)	0.11 / 4.33							EN455-2
Minimum Single Wall Finger Thickness (mm/mil)	0.14 / 5.51							EN 420:2003+A1:2009 / ASTM D3767
Minimum Single Wall Cuff Thickness (mm/mil)	0.09 / 3.54							EN 420:2003+A1:2009 / ASTM D3767
	During Ageing							
Force at Break (N)	≥ 9 N							EN 455-2

IONIC CONTENT

Concentration in µg/cm ²	Typical	Upper	Concentration in µg/cm ²	Typical	Upper
Ammonium	0.005	<0.06	Nitrate	0.795	<1.0000
Bromide	Not Detected	<0.0004	Nitrite	Not Detected	<0.0004
Calcium	0.582	<1.250	Phosphate	Not Detected	<0.0004
Chloride	0.494		Potassium	0.378	<0.080
Fluoride	Not Detected	<0.0004	Sodium	0.123	<0.150
Lithium	Not Detected	<0.0004	Sulphate	0.041	<0.150
Magnesium	0.025	<0.050	Zinc	Not Tested	<0.030

ORDERING INFORMATION

	SIZE	6	6.5	7	7.5	8	8.5	9
BUPS	REORDER NO.	BUPS-60	BUPS-65	BUPS-70	BUPS-75	BUPS-80	BUPS-85	BUPS-90

Technology



Performance Standards and Regulatory Compliance



For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV
Riverside Business Park
Blvd International, 55
1070 Brussels, Belgium
T: +32 (0) 2 528 74 00
F: +32 (0) 2 528 74 01

Asia Pacific Region

Ansell Global Trading Center
(Malaysia) Sdn Bhd
Prima 6, Prima Avenue
Block 3512, Jalan Teknokrat 6
63000 Cyberjaya
Selangor, Malaysia
T: +603 8310 6688
F: +603 8310 6699

North America Region

Ansell Healthcare Products LLC
111 Wood Avenue South,
Suite 210
Iselin, NJ 08830, USA
T: +1 800 800 0444
F: +1 800 800 0445

Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.
Bldv. Bernardo Quintana No. 7001-C,
Q7001 Torre II,
Suites 1304, 1305 y 1306,
Col. Centro Sur, c.p. 76079
Queretaro, Qro. Mexico
T: +52 442 248 1544 / 248 3133

Australia

Ansell Limited
Level 3, 678 Victoria Street,
Richmond, Vic, 3121
Australia
T: +61 1800 337 041
F: +61 1800 803 578

Russia

Ahcenn PYC
123610 Poccnr, Mockba
Kpachonpechehckar
Ha6. 12, n.3, op. 1304A
Ten. +7 495 258 13 16

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking © 2019 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

