

BIOCLEAN EXCELL

Sterile Nitrile Gloves

BEXS

About

BioClean Excell Sterile Nitrile Cleanroom Gloves are latex-free, resistance to a range of chemicals* and antistatic. Hand specific to reduce wearer hand fatigue, the gloves feature a textured surface enabling good grip and dexterity.

Specifications

COMPATIBILITY:

ISO Class 4 & EU GMP Grade A

LENGTH:

300mm (12")

MATERIAL:

Nitrile

PROTEIN LEVEL:

Nitrile contains no natural latex proteins

SURFACE:

Textured

SHAPE:

Hand specific

COLOUR:

White

STERILIZATION:

Gamma irradiation

CHEMICAL RESISTANT*

ACCELERATOR-FREE

ANTISTATIC

Features

- Resistant to a range of chemicals*
- Accelerator-free
- Powder-free & latex-free
- Good ESD properties
- Non-particulating EasyTear packaging



Ms.G.Love© represents BioClean's range of cleanroom gloves



Sterilization

- Method: Gamma irradiation
Minimum Dose: 25kGy
Sterility Assurance Level: 10^{-6}

Five (5) years from date of manufacture.
Store in a dry, cool place (<40°C) away
from direct sunlight and fluorescent
light.

CHARACTERISTICS	VALUE	TEST METHOD
Freedom From Holes	1.5 AQL Performance Level 2	EN 374-2:2003
Glove Length +/-10mm (mm/inches)	290/11.4	EN 420:2003
Single Wall Thickness Minimum (mils/mm)	Cuff: 3.54/0.09 Palm: 4.72/0.12 Finger Tip: 6.69/0.17	EN 455-2:2009
Force At Break During Shelf Life	≥ 9 N	EN 455-2:2009

PALM WIDTHS (mm) $\pm$ 5mm								
SIZE	6.0	6.5	7.0	7.5	8.0	8.5	9.0	10.0
mm	77	83	89	95	102	108	114	121
TEST METHOD					EN 420:2003			

PARTICLES					TYPICAL PARTICLE COUNT			TEST METHOD
≥ 0.5µm(counts/cm ²)					3000			IENT-RP-CC005.4
TYPICAL EXTRACTABLE IONS								TEST METHOD
ANIONS	Fluoride	Chloride	Nitrite	Bromide	Nitrate	Phosphate	Sulphate	IENT-RP-CC005.4
(µg/cm ²)	0.001	0.168	NT	NT	1.979	ND	0.325	
CATIONS	Lithium	Sodium	Ammonium	Potassium	Calcium	Magnesium	Zinc	
(µg/cm ²)	NT	0.314	NT	0.165	1.250	ND	0.050	
ND = Not Detected, NT = Not Tested, NA = Not Applicable								

TO ORDER

Re-order code	Size	Packing
BEXS-60	6.0	One pair per inner PE wallet; one wallet per sealed EasyTear PE pouch; 10 pouches per sealed outer PE bag; 20 outer bags per lined carton (200 pairs)
BEXS-65	6.5	
BEXS-70	7.0	
BEXS-75	7.5	
BEXS-80	8.0	
BEXS-85	8.5	
BEXS-90	9.0	
BEXS-100	10.0	

For additional information visit us at www.bioclean.com, or call us at +44 1638 663338

Nitritex Ltd

Minton Enterprise Park
Oaks Drive
Newmarket, Suffolk
CB8 7YY
United Kingdom

BEXS/1806/PDS28



Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

*Warning: No glove provides complete protection against all chemicals. Users must test the gloves against the particular chemicals and environment in which they will be used.

**Please see product validation pack or contact Ansell customer service for specific data on use of gloves with cytotoxic drugs. Gloves used for protection against such drugs must be selected specifically for the type of chemicals used.

Ansell, ® and TM are trademarks owned by Ansell Limited or one of its affiliates. © 2012 – 2018 Ansell Limited. All Rights Reserved.