

Magic Touch™ Insulated Labware

Keep Your Cool with Efficiency and Style!

Lightest, most durable icewares on the market, Magic Touch™ Ice Buckets and Pans have excellent insulating properties for keeping temperature sensitive samples in an optimal environment. Circular and rectangular floating racks fit neatly inside buckets and pans.



Expanded Polymer Foam Performs!

- Non-sweating, superior insulation properties
- Temperature range: -196° to 93°C (-321° to 200°F)
- Lightest, most durable insulated labware on the market
- Impervious to moisture and odors
- Chemical resistant
- Unbreakable under normal use conditions
- Use with ice, dry ice, dry ice alcohol slurries or liquid nitrogen
- Use multiple colors to track and differentiate samples

Four vibrant mix and match colors!



* As a service to our international dealers 12-pack bulk case numbers have been created for Magic Touch™ ice buckets and lab pans. Bulk case pack items are individually bagged but are not separately boxed in order to save you on freight costs.

Ice Buckets with Lids

Lightweight Easy Transport



A. Handy, 2.5L



B. N'Icer, 4.0L

- Non-sweating, superior insulation properties
- Use with ice, dry ice, dry ice alcohol slurries or liquid nitrogen
- Unbreakable under normal use conditions
- Flared edge for easy carrying
- Lid can be cut to make openings for thermometers, cooling coils, flasks etc.

Two Sizes:

A. 2.5 Liter, Handy

22.2cm diameter rim x 12.0cmH (8¾ x 4¾")

B. 4.0 Liter, N'Icer

21.5cm diameter rim x 17.0cmH (8½ x 6¾")



EACH* CATALOG NO.	PACK OF 12* CATALOG NO.	DESCRIPTION
M18848-2001	M18848-2091	Blue, 2.5L Handy with lid
M18848-2002	M18848-2092	Black, 2.5L Handy with lid
M18848-2003	M18848-2093	Red, 2.5L Handy with lid
M18848-2004	M18848-2094	Green, 2.5L Handy with lid
M18848-4001	M18848-4091	Blue, 4.0L N'Icer with lid
M18848-4002	M18848-4092	Black, 4.0L N'Icer with lid
M18848-4003	M18848-4093	Red, 4.0L N'Icer with lid
M18848-4004	M18848-4094	Green, 4.0L N'Icer with lid

Magic Touch™ Round Floating Tube Racks, See Page 187

* See page 184 for information regarding special case packing for international shipments

Insulated Lab Pans

Lightweight with Largest Surface Area

- Use with ice, dry ice, dry ice alcohol slurries or liquid nitrogen
- Lightweight with superior insulating properties
- Stackable, non-sweating and leakproof



A. Mini



B. Midi



C. Maxi



Three Sizes:

A. 1.0 Liter, Mini - 15.2 x 15.2 x 10.2cmH (6 x 6 x 4")

B. 4.0 Liter, Midi - 33 x 25.4 x 12.7cmH (13 x 10 x 5")

Holds 1 rectangular floating tube rack (see opposite page)

C. 9.0 Liter, Maxi - 40.6 x 33 x 12.7cmH (16 x 13 x 5")

Holds up to 3 rectangular floating tube racks (see opposite page)

EACH* CATALOG NO.	PACK OF 12* CATALOG NO.	DESCRIPTION
M18848-1101	M18848-1191	Blue, 1.0L Mini
M18848-1102	M18848-1192	Black, 1.0L Mini
M18848-1103	M18848-1193	Red, 1.0L Mini
M18848-1104	M18848-1194	Green, 1.0L Mini
M18848-4101	M18848-4191	Blue, 4.0L Midi
M18848-4102	M18848-4192	Black, 4.0L Midi
M18848-4103	M18848-4193	Red, 4.0L Midi
M18848-4104	M18848-4194	Green, 4.0L Midi
M18848-9101	M18848-9191	Blue, 9.0L Maxi
M18848-9102	M18848-9192	Black, 9.0L Maxi
M18848-9103	M18848-9193	Red, 9.0L Maxi
M18848-9104	M18848-9194	Green, 9.0L Maxi

* See page 184 for information regarding special case packing for international shipments

Floating Tube Racks

Convenient to Use, Easy to Store and Transport

- Tubes go in easy, expanded polymer foam holds them secure, yet they come out easily when needed
- Racks float in water, dry ice alcohol slurries or liquid nitrogen
- Use in buckets and pans





Round Floating Tube Racks

Fit N'icer and Handy Ice Buckets

21cm diameter (8¼") • Two Hole Sizes

- A. 1.5ml microcentrifuge tubes,**
44 places – 1cm thick (⅜")
- B. 15ml conical tubes,**
30 places – 1.58cm thick (⅝")



Rectangular Floating Tube Racks

Fit Midi and Maxi Ice Pans

23.5 x 10.8cm (9¼ x 4¼") • Two Hole Sizes

- C. 1.5ml microcentrifuge tubes,**
36 places – 1cm thick (⅜")
- D. 15ml conical tubes,**
24 places – 1.58cm thick (⅝")



Order by color suffix in second table for complete catalog numbers

Per each, 36 per case.

CATALOG NO	DESCRIPTION	PLACES
A. M18848 -	Round Rack for 1.5ml Tubes	44
B. M18848 -	Round Rack for 15ml Tubes	30
C. M18848 -	Rectangular Rack for 1.5ml Tubes	36
D. M18848 -	Rectangular Rack for 15ml Tubes	30

Color suffix listing

CATALOG NO	BLUE	BLACK	RED	GREEN
A. M18848 -	0011	0012	0013	0014
B. M18848 -	0021	0022	0023	0024
C. M18848 -	0111	0112	0113	0114
D. M18848 -	0121	0122	0123	0124