



Write-On™ Multi-Roll Label Tape Dispensers

Mix Tape Colors and Widths to Suit Your Labeling Needs

U.S. Patent 6,974,060

Write-On™ Label Tape Dispensers provide a smooth writing surface, making labeling easy and fast. When marked tape is advanced and cut off on the serrated edge, a fresh length of label tape is positioned on the writing platform.

Various width and color tape rolls can be mixed together on the spindles, which easily accommodate 5 or more 12 - 19mm ($\frac{1}{2}$ - $\frac{3}{4}$ ") narrow rolls or 4 - 25mm (1") wide rolls. The spindles are easily removed from the unit for loading and unloading, and can be adjusted to hold either 25 or 75mm (1 or 3") core diameter rolls. Dispenser includes three plastic separating discs to assure smooth, snag free dispensing for up to four rolls of tape. Additional separating discs are available separately (see accessories). Molded-in pen holder wells on the dispenser help to keep writing instruments handy. Ideal for use with the full line of F13463 series Write-On™ Tapes, which are available in 6 colors and three widths; 12, 19, and 25mm ($\frac{1}{2}$, $\frac{3}{4}$ and 1").

Two models are available:

- A. Bench-Top / Desk Dispenser** has a plastic coated, weighted steel base with rubber feet to prevent slipping on the benchtop. The serrated cutting edge is steel for long term use. Includes one free roll of $\frac{1}{2}$ ", White, Write-On™ Labeling Tape – a \$9.00 value!
- B. Economy Dispenser** is designed for lighter duty use as the plastic cutting edge is molded into the housing. Includes mounting holes for secure attachment to a benchtop (screws not supplied). Can also be attached to vertical surfaces using Wall Mount adapter-F13473-0001 (opposite page). Includes one free roll of $\frac{1}{2}$ " White, Write-On™ Labeling Tape – a \$9.00 value!

CATALOG NO.	DESCRIPTION	DIMENSIONS (L x W x H)
F13472-0000	A. Benchtop Model	30.2 x 14.3 x 12.1cm (12 x 5 $\frac{5}{8}$ x 4 $\frac{3}{4}$ ")
F13473-0000	B. Economy Model	30.2 x 14.3 x 10.5cm (12 x 5 $\frac{5}{8}$ x 4 $\frac{1}{4}$ ")

Wall Mount Adapter for Economy Write-On™ Label Tape Dispenser

Economy Model Dispenser (F13473-0000) can be wall mounted vertically or horizontally with the platform positioned to the left or right. The U-shaped stainless steel adapter presents the dispenser at a 45° angle off the wall to facilitate writing. Wet Surface Pen (page 193) offers consistent ink flow on vertical surfaces.



CATALOG NO.	DESCRIPTION	DIMENSIONS
F13473-0001	Wall Mount Adapter	11 x 13cm (4½ x 5½")

Separating Discs for Write-On™ Label Tape Dispensers

Separating discs improve the function of the dispensers when loaded with more than four rolls of labeling tape. **3 per bag.**

CATALOG NO.	DESCRIPTION
F13469-0003	Separator Discs



Write-On™ Single Roll Label Tape Dispenser

Make Your Own Stick-On Labels

Convenient dispenser holds and releases tape cut to your specification. Writing platform accommodates both left and right handed users for labeling tape prior to dispensing. The epoxy-coated steel platform has finger holes for releasing the tape and a cut off edge for one hand use. A weighted base with non-skid pads prevents sliding or lifting. The 7.6cm (3") core fits standard tape rolls up to 25mm (1") wide and accepts tape rolls having a maximum diameter of 13.3cm (5½"). The base is 21.3 x 8.7 x 9.5cm H (8½ x 3½ x 3¾") with a 10.2 x 15.2cm (4 x 6") platform. Shipping weight 1.8kg (4lb). Accepts F13463 series tapes in both 40 yard and 60 yard rolls. Tape not included.

CATALOG NO.
H13461-0000



Write-On™ Label Tape

Available in 6 Vivid Colors. Great for Use in Low Temperature Applications

Adheres to any clean surface, including Teflon® fluoropolymer resin coated materials, and peels off without leaving a sticky residue. Pencil, ball-point pen or solvent ink marker remain legible under the most demanding laboratory conditions. Consists of colored paper with rubber resin adhesive that resists moisture and most solvents. The tape will withstand a maximum temperature of 157°C (315°F) for 30 minutes or 80°C (176°F) for 24 hours and can be used down to -80°C (-112°F). Withstands a steam autoclaving cycle at 121°C (250°F) and gamma irradiation between 13.0 and 45.0kGy. Excellent for use with Write-On™ Multi-Roll Label Tape Dispensers. **Per each, 24 rolls per case.**



CATALOG NO.	WIDTH	CORE	LENGTH	COLOR
F13462-0005	12.7mm (½")	2.54cm (1")	13.7m (15yds)	white
F13462-0075	19.0mm (¾")	2.54cm (1")	13.7m (15yds)	white
F13462-0010	25.4mm (1")	2.54cm (1")	13.7m (15yds)	white
F13463-0005	12.7mm (½")	7.6cm (3")	36.4m (40yds)	white
F13463-0075	19.0mm (¾")	7.6cm (3")	36.4m (40yds)	white
F13463-0010	25.4mm (1")	7.6cm (3")	36.4m (40yds)	white
F13463-0105	12.7mm (½")	7.6cm (3")	54.6m (60yds)	white
F13463-0175	19.0mm (¾")	7.6cm (3")	54.6m (60yds)	white
F13463-0110	25.4mm (1")	7.6cm (3")	54.6m (60yds)	white
F13463-1005	12.7mm (½")	7.6cm (3")	36.4m (40yds)	red
F13463-1075	19.0mm (¾")	7.6cm (3")	36.4m (40yds)	red
F13463-1010	25.4mm (1")	7.6cm (3")	36.4m (40yds)	red
F13462-2005	12.7mm (½")	2.54cm (1")	13.7m (15yds)	yellow
F13462-2075	19.0mm (¾")	2.54cm (1")	13.7m (15yds)	yellow
F13462-2010	25.4mm (1")	2.54cm (1")	13.7m (15yds)	yellow
F13463-2005	12.7mm (½")	7.6cm (3")	36.4m (40yds)	yellow
F13463-2075	19.0mm (¾")	7.6cm (3")	36.4m (40yds)	yellow
F13463-2010	25.4mm (1")	7.6cm (3")	36.4m (40yds)	yellow
F13462-3005	12.7mm (½")	2.54cm (1")	13.7m (15yds)	green
F13462-3075	19.0mm (¾")	2.54cm (1")	13.7m (15yds)	green
F13462-3010	25.4mm (1")	2.54cm (1")	13.7m (15yds)	green
F13463-3005	12.7mm (½")	7.6cm (3")	36.4m (40yds)	green
F13463-3075	19.0mm (¾")	7.6cm (3")	36.4m (40yds)	green
F13463-3010	25.4mm (1")	7.6cm (3")	36.4m (40yds)	green
F13463-4005	12.7mm (½")	7.6cm (3")	36.4m (40yds)	blue
F13463-4075	19.0mm (¾")	7.6cm (3")	36.4m (40yds)	blue
F13463-4010	25.4mm (1")	7.6cm (3")	36.4m (40yds)	blue
F13463-5005	12.7mm (½")	7.6cm (3")	36.4m (40yds)	orange
F13463-5075	19.0mm (¾")	7.6cm (3")	36.4m (40yds)	orange
F13463-5010	25.4mm (1")	7.6cm (3")	36.4m (40yds)	orange

Write-On™ Label Tape Multi-Pack

3-pack of tape in Green, Yellow and White. Each roll is 1.9cm ($\frac{3}{4}$ ") wide and 13.7m (15 yards) long with a 2.54cm (1") inner core. Use Green for safety labels, Yellow for chemical reactivity labels and White for specific hazards/labels.



CATALOG NO.

F13462-0300

DESCRIPTION

1 Each Green, Yellow, White
 $\frac{3}{4}$ " W x 15 yards, 1" core



Write-On™ Label Tape Rainbow 6-Pack

Six tape colors that correspond with familiar hazard color codes to highlight safety information written on the labels. Use Yellow for chemical reactivity, Green for safety labels, Blue indicates health hazards, Red alerts to flammability; Orange indicates a biohazard and White can describe a specific hazard. Each roll is 1.9cm ($\frac{3}{4}$ ") wide and 36.4cm (40 yards) long with a 7.6cm (3") inner core.

CATALOG NO.

F13463-0600

DESCRIPTION

1 Each Green, White, Yellow, Red, Orange, Blue
 $\frac{3}{4}$ " W x 40 yards, 3" core

Multi-Roll Tape Dispenser

Flat Writing Surface Helps Make Neater Labels

Scienceware® Write-On™ Labeling Tape or any tape with a 7.6cm (3") ID core will fit on the drum of this multi-roll dispenser, which includes disks to prevent rolls from binding against each other. Large writing surface, 10.2 x 15.2cm (4" x 6") for preparing tape labels, has a built in serrated cutting edge. **Per each.**



CATALOG NO.

H13468-0000

DIMENSIONS (L x W x H)

26 x 15 x 11cm (10 $\frac{1}{4}$ x 6 x 4 $\frac{3}{8}$ ")