


## BEAKERS, LOW FORM, WITH SPOUT

Kádinka nízká s výlevkou

Becher, niedrige Form, mit Ausguß

Béchers, forme basse à bec

Vasos de precipitación graduados, con pico, forma baja



Cat. No. 155	ISO 3819, DIN 12 331, CSN 70 4031					
	Code	Capacity	d	h	Packing	EAN
	600000632	ml	mm	mm	Quantity	Carton
155/5	411010005	5	22	30	10	8593419605819
155/10	411010010	10	26	35	10	8593419605833
155/25	411010025	25	34	50	10	8593419602436
155/50	411010050	50	42	60	10	8593419639036
155/100	417010100	100	50	70	10	8593419602467
155/150	417010150	150	60	80	10	8593419602474
155/250	417010250	250	70	95	10	8593419602481
155/400	417010400	400	80	110	10	8593419602498
155/600	417010600	600	90	125	10	8593419602504
155/800	417010800	800	100	135	10	8593419602511
155/1000	417010940	1000	105	145	10	8593419602528
155/2000	417010950	2000	130	185	4	8593419602535
155/3000	411010952	3000	150	210	4	8593419602542
155/5000	411010956	5000	170	270	2	8593419602559
155/10000	411010966	10000*	217	350	1	8593419646119


## BEAKERS, TALL FORM, WITH SPOUT

Kádinka vysoká s výlevkou

Becher, hohe Form, mit Ausguß

Béchers, forme haute à bec

Vasos de precipitación graduados, con pico, forma alta



Cat. No. 153	ISO 3819, DIN 12 331, CSN 70 4031					
	Code	Capacity	d	h	Packing	EAN
	600000632	ml	mm	mm	Quantity	Carton
153/25	411012025	25	30	55	10	8593419605857
153/50	411012050	50	38	70	10	8593419639081
153/100	411012100	100	48	80	10	8593419602573
153/150	411012150	150	54	95	10	8593419602580
153/250	417012250	250	60	120	10	8593419602597
153/400	417012400	400	70	130	10	8593419602603
153/600	417012600	600	80	150	10	8593419602610
153/800	417012800	800	90	175	10	8593419602627
153/1000	417012940	1000	95	180	10	8593419602634
153/2000	417012950	2000	120	240	6	8593419602641
153/3000	411012952	3000	135	280	6	8593419602658

\* This size is out of norm ISO, DIN. / Tato velikost je mimo normu ISO, DIN. / Sondergrosse, nicht nach ISO, DIN  
 Taille spéciale, non conforme aux normes ISO, DIN / Esta dimensión no corresponde a las normas ISO, DIN