

SKH-MINI... **A fully disposable spill kit to clean up small chemical spills quickly, easily and without risk of injuries.**

- Highly visible green colour.
- The autoclavable disposal bag is heat resistant up to 145°C and therefore guarantees a safe disposal of laboratory waste.
- Adsorbs chemicals, dangerous liquids, water & oil based fluids.
- Suited for spills up to 1.7 litres.
- The specially designed disposable brush is intended to sweep up any broken glass or debris allowing a safe clean up without risk of injuries.
- Thanks to the compact design, the SKH-MINI can easily be placed in areas where spills are most likely to happen.
- Protective gloves are powder free, have roughened fingertips, increased tear resistance and a good chemical resistance for more safety.



Order Ref.	Contents	Packaging	Adsorption Capacity
SKH-MINI	2 pads 41 cm x 51 cm 1 pair of disposable chemical protective gloves 1 autoclavable disposal bag 1 disposable brush 1 multilingual instruction guide	4 kits/box	1.7 litre/kit