


## Sampling bag SteriBag | the sterilised packaging bag



The sterilised sampling bag SteriBag for samples in chemistry, pharmaceuticals, medicine, veterinary science and foodstuffs etc.


The bags are made of polyethylene and are therefore exceptionally tearproof and leakproof.

The SteriBag sampling bag is distinguished by its safe handling. When used according to instructions, samples are guaranteed sterile.

Integrated safety: The tab-ends are sealed in plastic film to prevent skin injuries.

- ▶ Sterilised
- ▶ Can be labelled
- ▶ Leakproof
- ▶ Save space
- ▶ Transparent
- ▶ With practical closure to tear-open
- ▶ Proven millions of times

**Bulk prices**



Sampling bag SteriBag											
Capacity ml	Model	L x W mm	Pieces per pack	Price €	Quantity scale from	Scale price €	Item no.				
60	without labelling area	127 x 76	500	112.50	6	106.55	5344-0060	N3			
60	with labelling area	127 x 76	500	122.25	6	116.20	5344-0061	N3			
150	without labelling area	178 x 76	500	116.20	6	108.90	5344-0180	N3			
150	with labelling area	178 x 76	500	125.90	6	119.75	5344-0181	N3			
450	without labelling area	229 x 114	500	136.70	6	130.70	5344-0540	N3			
450	with labelling area	229 x 114	500	145.15	6	136.70	5344-0541	N3			
650	without labelling area	229 x 140	500	152.45	6	143.95	5344-0720	N3			
650	with labelling area	229 x 140	500	174.20	6	166.95	5344-0721	N3			
900	without labelling area	305 x 114	500	158.55	6	151.25	5344-0810	N3			
900	with labelling area	305 x 114	500	168.15	6	159.75	5344-0811	N3			
1650	without labelling area	305 x 178	250	121.00	12	116.20	5344-1650	N3			
1650	with labelling area	305 x 178	250	136.70	12	125.90	5344-1651	N3			

## SteriBag StandUp sample bag

SteriBag StandUp is a sterile sample bag with a stand-up base. The special bottom unfolds as soon as the bag is filled. This makes SteriBag StandUp stable and able to stand on its own without a separate stand or support, leaving both your hands free to take the sample and fill the bag. The wire bracket at the closure keeps the bag open, making it easy and convenient to fill.

The sample bag is the ideal space-saving, unbreakable, inexpensive and sterile replacement for bottles and beakers.

Can be used for samples available as liquids, powders, solids and pastes; suitable for foods and can be frozen.

The stable standing SteriBag StandUp sample bag is easy and safe to transport.

SteriBag StandUp, the sample bag for easy handling of samples.

- ▶ Stands stably
- ▶ Protected against discharge
- ▶ Sterilised
- ▶ With labelling field



**Bulk prices**

### SteriBag StandUp sample bag

Capacity ml	L x W mm	Pieces per pack	Price €	Quantity scale from	Scale price €	Item no.
400	229x114	500	179.25	5	170.25	5344-0400 N3
600	229x152	500	187.35	5	177.95	5344-0600 N3