



TH 15 mit KS 15 und Aufsatz Combifix KS
TH 15 with KS 15 and rack system Combifix KS

TH 15 Nr. / No. 6161 000

Beschreibung

Konstante Temperaturen; auch für hohe Gefäße
(z.B. 1000 ml Erlenmeyerkolben)
Passend auf alle Kompaktschüttler oder als eigenständiges Gerät

Lieferumfang Grundgerät

Beschichtetes Metallgehäuse mit 3 Fenstern, Fronttür nach oben zu öffnen

Technische Daten

Temperaturbereich: 5°C über Raumtemperatur
bis +50°C
Temperaturkonstanz: +/- 1% vom Einstellwert
Innenabmessungen (B x T x H): 490 x 410 x 310 mm
(Höhe ausreichend für 1000 ml Erlenmeyerkolben)
Temperaturfühler: PT 100
Luftumwälzung: 240 m³/h
Elektrischer Anschluss: 230 V, 50 / 60 Hz
Ausführung 115 V auf Anfrage
Außenmaße (B x T x H): 510 x 450 x 455 mm
Gewicht: 21 kg

Sicherheitsabschaltung des Schüttlers in Kombination mit Kompaktschüttlern
Standard-Regler 5° über RT bis +50°C

Option

Erweiterter Temperaturbereich bis +60°C

Nr. 0052 105

Einbau einer Zusatzheizung, um Temperaturen bis +60°C (max.)
zu erreichen. Einbau bzw. Nachrüstung nur im Werk möglich.

Description

Constant temperatures; also for high vessels
(e.g. 1000 ml Erlenmeyer flasks)
Can be combined with all compact shakers or used as independent unit

Basic equipment

Coated metal housing with 3 windows, the front door opens upwards

Technical Data

Temperature range: 5°C above ambient
to +50°C
Temperature accuracy: +/- 1% of the set value
Inner dimensions (w x d x h): 490 x 410 x 310 mm
(height sufficient for 1000 ml Erlenmeyer flasks)
Temperature sensor: PT 100
Air circulation: 240 m³/h
Electrical supply: 230 V, 50 / 60 Hz
115 V version on request
External Dimensions (w x d x h): 510 x 450 x 455 mm
Weight: 21 kg

Safety cut-off of the shaker if combined with a Bühler Compact Shaker
Standard temperature controller 5°C above ambient to +50°C

Option

Extended temperature range up to +60°C

No. 0052 105

Installation of an additional heater in order to reach temperatures up
to max. +60°C. Installation only possible in our works.



Einhängetablar TH 15
Nr. 0052 126

Zusätzliche Ablage für Proben,
die nur temperiert werden müssen.
Abmessungen: 378 x 318 mm.

Shelf TH 15
No. 0052 126

Additional tray for samples which
have to be temperature controlled
but need no shaking.
Dimensions: 378 x 318 mm.