

LABORHITS

Detectable disposable sampler, SteriPlast®



- Due to special additive in the material, they can be detected by colour scanners, metal detectors or X-rays
- PS detectable, blue
- Individually packaged and sterilised by gamma rays
- Production, assembly and packaging according to clean room class 7
- Complies with EU food requirements and FDA regulations




Detectable scoops, blue

The food scoop is ideal for powders, granulates and pastes.

Capacity ml	Overall length mm	PK	Cat. No.
25	141	10	6.284 563
50	160	10	6.284 564
100	205	10	6.284 565
150	216	10	6.284 566
250	232	10	6.284 567
500	280	10	6.284 568
1000*	332	10	6.284 569

* not produced in class 7 clean room



Detectable spatulas, blue

As a result of the long, ergonomically-shaped, stable handle and a sharp blade, it is possible to even penetrate directly into containers such as paper or plastic sacks. The spatula for foodstuffs is ideal for sampling powders, granulates and pastes.

Lenght mm	Imm. depth mm	PK	Cat. No.
263	150	10	6.284 570



Detectable spoon, blue

The spoon for foodstuffs is available in two sizes, 2.5 ml (similar to a teaspoon) and 10 ml (similar to a tablespoon), and is ideal for sampling powders, granulates, pastes and fluids.

Capacity ml	Lenght mm	PK	Cat. No.
2.5	127	100	6.284 571
10.0	170	100	6.284 572



Detectable scraper, blue

Sturdy and convenient scraper for disposable use. The scraper is suitable for scraping off, smoothing or cleaning as well as removing residue.

Width mm	Length mm	PK	Cat. No.
80	200	10	6.284 573