

Customer information

05.06.2023, Version 2, FRI

NANOCOLOR® VIS II – Analysis without limits

Flexible power supply via power cable, battery and car's electrical system

Measurements are no longer a barrier on-site with the spectrophotometer NANOCOLOR® VIS II due to the flexible power supply. Beside the power supply via the supplied power cable the device can also be powered via the voltage of the car's electrical system together with the car's adapter cable or an external battery. We have tested its functionality in combination with the battery model "Dino start-up support 136102". The compatibility is confirmed for the following specifications:

Specifications

Description	Dino start-up support Art.-No:136102
Battery charge	Fully charged before use
Temperature range	20 °C – 45 °C



NANOCOLOR® VIS II



Dino start-up 136102



Car adapter cable

Description

Spectrophotometer NANOCOLOR® VIS II
incl. manual (quick start guide), touchpen, protective covering, power cable, USB cable, USB-stick,
calibration cuvette, certificate and cleaning cloth in cardboard box

REF

919650

Car adapter cable for NANOCOLOR® VIS II

919938

External battery for NANOCOLOR® VIS II

Available online via different suppliers

Contact

If you have further questions, please do not hesitate to contact us:

Technical support and customer service

Tel.: +49 24 21 969-333

E-Mail: support@mn-net.com

For orders (international)

E-Mail: sales@mn-net.com


Management
System
EN ISO 13485:2016
ISO 9001:2015

www.tuv.com
ID 0000056401

MACHEREY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Germany · www.mn-net.com

DE +49 24 21 969-0 info@mn-net.com

CH +41 62 388 55 00 sales-ch@mn-net.com

FR +33 388 68 22 68 sales-fr@mn-net.com

US +1 888 321 62 24 sales-us@mn-net.com