



Certificate of Accuracy and Reproducibility

LLG - Microlitre Pipettes with fixed volume				
Capacity μl	Accuracy $\leq \pm A\%$ (max. vol.)	Reproducibility $\leq \pm CV\%$ (max. vol.)	Colour code	Cat. No.
5	1.0	0.5	grey	9.277 515
10	1.0	0.5	grey	9.277 516
20	4.5	1.0	grey	9.277 517
25	2.5	0.6	yellow	9.277 518
50	2.5	0.6	yellow	9.277 519
100	2.5	0.6	yellow	9.277 520
200	2.5	0.6	blue	9.277 521
250	2.9	0.6	blue	9.277 522
500	2.9	0.6	blue	9.277 523
1000	2.4	0.6	violet	9.277 524
2000	2.4	0.6	violet	9.277 525
2500	2.4	0.6	violet	9.277 526
5000	2.4	0.6	turquoise	9.277 527
10.000	3.0	0.6	turquoise	9.277 528

LLG - Microlitre Pipettes with variable volume				
Capacity μl	Accuracy $\leq \pm A\%$ (max. vol.)	Reproducibility $\leq \pm CV\%$ (max. vol.)	Colour code	Cat. No.
0.1 to 2,5	1.5	0.8	grey	9.280 000
0.5 to 10	1.0	0.5	grey	9.280 001
2 to 20	4.5	1.0	grey	9.280 002
10 to 100	2.8	0.8	yellow	9.280 003
20 to 200	2.5	0.6	yellow	9.280 004
100 to 1000	2.9	0.6	blue	9.280 005
500 to 5000	2.4	0.6	purple	9.280 006
1000 to 10.000	3.0	0.6	turquoise	9.280 007

Lab Logistics Group GmbH
Product Management
Am Hambuch 1
D-53340 Meckenheim
Germany

Meckenheim, 2014-10-13

This document was sent online and is valid without signature

