


Comparison table Flocculators



		JLT 4	JLT 6	FC4S	FC6S	FP4
Number of positions		4	6	4	6	4 (portable)
Interface		Digital	Digital	Analog	Analog	Digital
Analysis condition		Same speed for each rod	Same speed for each rod	Independent for each rod	Independent for each rod	Same speed for each rod
Speed settings		From 10 to 300 rpm	From 10 to 300 rpm	10-15-30-45-60-90-120-150-200-300	10-15-30-45-60-90-120-150-200-300	From 10 to 300 rpm
Time settings		1÷999 min or 1÷99 hours Continuous	1÷999 min or 1÷99 hours Continuous	- Continuous	- Continuous	1 min ÷ 99h 59 min Continuous
Lighting system		Illuminated back panel	Illuminated back panel	Illuminated back panel	Illuminated back panel	Central light
Power input		100÷240 V	100÷240 V	100÷240 V	100÷240 V	100÷240 V
Power		19 W	19 W	26 W	32 W	6 W
Weight	kg	13.0	17.0	12.5	18.0	3.8
	lb	(28.6)	(37.4)	(27.5)	(39.6)	(8.3 lb)
Dimensions (WxHxD)	mm	655x404x296	935x404x296	655x404x296	935x404x296	250x310x265
	in	(25.7x15.9x11.6)	15.2x30.7x16.4	(25.7x15.9x11.6)	15.2x30.7x16.4	(9.8x12.2x10.4)