



### Corning® 15mL PET Centrifuge Tubes, Rack Packed with Plug Seal Cap, Sterile, 50/Rack, 500/Case (Product #430055)

Corning® 15mL clear polyethylene terephthalate (PET) centrifuge tubes have traditional polyethylene plug seal caps and are packed 50 per rack. Please check our Technical Information Listing for information on the physical and chemical properties of Corning plastic centrifuge tubes as well as important suggestions for safe use.

- Black printed graduations and a large white marking spot
- Maximum RCF is 3,600xg
- Sterile, certified nonpyrogenic and DNase-/RNase-free
- 95 k PA (14psi) pressure tested

### Product Details At-a-Glance

Qty./Pk.:	50
Qty./CS.:	500
Capacity:	15 mL
Tube Style:	Centrifuge
Bottom Style:	Conical
Neck Style:	Wide
Top Style:	Threaded
Max RCF:	3600 xg
Graduations:	Yes
Graduation Range:	1.5-14.5 mL
Graduation Interval:	0.5 mL
Color:	Clear
Cap Style:	Plug Seal
Cap Material:	High Density Polyethylene
Cap Diameter:	22.9 mm (approx)
Cap Color:	Orange
Sterile:	Yes
Self-Standing:	No
Outer Diameter:	17.5 mm (approx)
Length:	118.5 mm (approx)
Feature:	Disposable
Rack:	Yes