



Corning® 250mL PP Centrifuge Tubes with Plug Seal Cap, Sterile, 6/Pack, 102/Case (Product #430776)

Corning® 250mL polypropylene (PP) centrifuge tubes are ideal for applications requiring large volume centrifugation. Please check our Technical Information Listing for information on the physical and chemical properties of Corning plastic centrifuge tubes as well as important suggestions for safe use.

- Support cushions (Catalog No. 430236 for 250) must be used with this product unless the centrifuge rotor has appropriate V-bottom cushions.
- Polyethylene plug seal caps.
- Maximum RCF of 6,000xg
- Tubes are sterile and certified nonpyrogenic.

Product Details At-a-Glance

Qty./Pk.:	6
Qty./CS.:	102
Capacity:	250 mL
Tube Style:	Centrifuge
Bottom Style:	Conical
Neck Style:	Short
Top Style:	Threaded
Max RCF:	6000 xg
Graduations:	Yes
Graduation Range:	1.5-14.5 mL
Graduation Interval:	0.5 mL
Color:	Natural
Cap Style:	Plug Seal
Cap Material:	High Density Polyethylene
Cap Color:	Orange
Sterile:	Yes
Self-Standing:	No
Feature:	Disposable
Rack:	No