

REF 01130

- A proven method of positively confirming a completed sterilization cycle in both gravity discharge and vacuum autoclaves. White indicator lines change to black when proper levels of moisture and temperature have been achieved. The uniquely printed lines cannot be rubbed off and will not contaminate surgical instruments or dressings.
- This easy-to-handle tape adheres on contact and stays in place through live steam pressure.
- This tape is made of a saturated crepe paper printed with white indicator lines that turn to black once autoclave process is completed. It is coated with a dry natural rubber based adhesive.
- Indications for Use: The Process Indicator Tape for Steam Sterilization is indicated for use in holding sterilization packs together and as process indicators in gravity steam sterilizers operating at 121°C and pre-vacuum steam sterilizers operating at 132-134°C. The white indicator stripes on the tape turn black when exposed to steam sterilization providing indication of processed items.

TECHNICAL DATA

	Natural with Printed White Indicator Lines
Backing	Saturated Crepe Paper
Adhesive	Dry Natural Rubber Based Adhesive
Thickness	6.0 mil (0.15mm)
Adhesion to steel	25oz/in (6.7N/25mm)
Tensile strength	20 lb/in (89N/25mm)
Elongation	6%
Standard Sterilization Cycle (with pre-vacuum cycle)	
	121 °C(250 °F) for 20 minutes
	132 to 134 °C (270 to 274°F) for 4 minutes
Meets Canadian Specification	CAN/CSA-Z11140-1-07 (R2012)
	ANSI/AAMI/ISO 11140-1:2005 (R2010)
	EN ISO 11140-1:2009

STOCK SIZE

18mm x 55m 48 rolls/case 24mm x 55m 36 rolls/case