

- Designed for processing blots and staining gels
- Optimal fixed speed and tilt
- Three dimensional motion
- Blot boxes perfect for protein gels



**Ordering Information**

Cat#	Description
S0600*	<b>Enduro MiniMix</b> 3-D shaker with 20 x 16.5 cm tray and non-slip rubber mat, 120V
E2110-BBS	Blotting box, 9.1 x 6.6 cm, 3-5 mL gel capacity, for mini protein gels, 10/pk
E2110-BBL	Blotting box, 11.7 x 8.9 cm, 6-10 mL gel capacity, for larger mini protein gels, 10/pk

\* To order 230V units add -230V to the end of the catalog number.  
To specify cord type add -EU or -UK

**Specifications:**

Speed range:	18 rpm (20 rpm for 230V) Fixed
Motion/Pitch:	3D/ Fixed 5°
Max load:	1.75lb/0.8 kg
Ambient oper range	+ 4° to 65° C
Dim. (WxDxH):	8 x 7 x 4.5 in 20.3 x 17.8 x 10.5 cm
Weight:	1.9 lb/0.88 kg



- Gyrotory action provides, gentle, thorough mixing
- Ideal for mixing samples in blood and centrifuge tubes
- Perfect for resuspension of large, fragile DNA fragments
- Included dimpled mat holds tubes in place
- Safe for incubator and cold room use

**Specifications:**

Speed range:	24 rpm (20 rpm for 230V)
Motion:	3D/ Fixed 20°
Max load:	1.75lb/0.8 kg
Ambient oper range	+ 4° to 65° C
Dim. (WxDxH):	9.5 x 6.75 x 6 in 24.1 x 17.15 x 15.24 cm
Weight:	1.9 lb/0.88 kg



**Ordering Information**

Cat#	Description
S0500*	<b>GyroMini</b> Nutating 3-D mixer with dimpled mat, 120 V
S0500-MAT	Dimpled mat (9.5 x 6.75)

\* To order 230V units add -230V to the end of the catalog number.  
To specify cord type add -EU or -UK