



designed for scientists



AS 501.3 Separating funnel attachment

/// Data Sheet

For shaking out, eluting, extracting, gassing out, dissolving, enriching, etc. Adjustment for the clamping rolls is infinitely variable, the set-up height can be changed by means of clamping devices.

Dimensions (W x D x H): 480 x 505 x 190 mm

Capacity: 4 x 500 or 3 x 1000 or 2 x 2000 ml

Weight: 3720 g



designed for scientists

Technical Data

Width [mm]	106
Depth [mm]	73.11
Weight [kg]	3