




designed for scientists



IKAFLON 15 round (5 pcs.)

/// Data Sheet

Dimensions (L x Ø): 15 x 6 mm

The round, cylindrical form without pivot ring ensures even mixing.

The PTFE coating is chemically inert and resistant to aggressive chemicals.

Suitable for a variety of vessels with flat bottoms.

Universal stirrer for a variety of applications.



designed for scientists

Technical Data

Length [mm]	15
Diameter [mm]	6
Shape	round
Coating	PTFE
Material in contact with medium	PTFE
Material: Magnet	Alnico V
Magnetic attraction force	normal
Working temperature max. [°C]	280
Weight [kg]	0.0091