




designed for scientists



IKAFLON 25 slide round (5 pcs.)

/// Data Sheet

Dimensions (L x Ø): 25 x 6 mm

Cylindrical magnetic stir bars provide outstanding centering and consistent mixing.

The integrated pivot ring reduces vibrations and wear during the stirring process.

The PTFE coating is chemically inert and resistant to aggressive chemicals

For use in containers with uneven or curved bottoms. The pivote ring maintains the optimum stirring position.



designed for scientists

Technical Data

Length [mm]	25
Diameter [mm]	6
Coating	PTFE
Material in contact with medium	PTFE
Material: Magnet	Alnico V
Magnetic attraction force	normal
Working temperature max. [°C]	280
Weight [kg]	0.015