




designed for scientists



IKAFLON 50 round (5 pcs.)

/// Data Sheet

Dimensions (L x Ø): 50 x 8 mm

The round, cylindrical form without pivot ring ensures even mixing.

The PTFE coating is chemically inert and resistant to aggressive chemicals.

Suitable for a variety of vessels with flat bottoms.

Universal stirrer for a variety of applications.



designed for scientists

Technical Data

Length [mm]	50
Diameter [mm]	8
Shape	round
Coating	PTFE
Material in contact with medium	PTFE
Material: Magnet	Alnico V
Magnetic attraction force	normal
Working temperature max. [°C]	280
Weight [kg]	0.047