



designed for scientists



LR 1000.65 pH Electrode receptacle

/// Data Sheet

Made of stainless steel with two FFKM seals to be connected to the port in the reactor cover instead of an NS 14/23 adapters using the M 24/1.5 thread.

Immersion length in the vessel: 13,5 cm

Total length of the adapter: 21 cm

With internal thread for connecting the pH sensor with the existing PG 13.5 adapter head.



designed for scientists

Technical Data

Dimensions (W x H x D) [mm]	430 x 175.5 x 55
Weight [kg]	0.3456