



designed for scientists



LR 2000.60 Sensor receptacle

/// Data Sheet

To accommodate the PT 100.30, PT 100.25, PT 1000.61 and H 68.55 temperature sensors with a diameter of 6 mm. Made of stainless steel with FFKM seals to be connected to the port in the reactor cover of the LR-2.ST instead of an NS 14/23 adapters using the M 24/1.5 thread.



designed for scientists

Technical Data

Dimensions (W x H x D) [mm]	1500 x 16 x 8
Weight [kg]	0.17