



designed for scientists



LR 2000.65 Receptacle for pH electrode

/// Data Sheet

Made of stainless steel with two FFKM seals to be connected to the port in the reactor cover instead of an NS 29/32 adapters using the M 36/1.5 thread.

Immersion length in the vessel: 14 cm

Total length of the adapter: 19 cm

With internal thread for connecting the pH sensor with the existing PG 13.5 adapter head.



designed for scientists

Technical Data

Depth [mm]	279.5
Height [mm]	4.74
Weight [kg]	0.4229