



MF 10.2 Impact grinding head

/// Data Sheet

For crushing brittle, hard materials such as minerals, building materials up to Mohs hardness 6. Before being discharged, the ground material passes through a sieve. This sieve is interchangeable and available in different hole sizes (not incl. with delivery). The ground material can then be collected using an NS 29 standard ground vessel.



designed for scientists

Technical Data

Process type	inline
Circumferential speed max. [m/s]	34
Feed hardness max. [Mohs]	6
Feed grain size max. [mm]	10
Material milling chamber	stainless steel 1.4308
Dimensions (W x H x D) [mm]	320 x 450 x 300
Weight [kg]	2.942
Permissible ambient temperature [°C]	5 - 40