



designed for scientists



MX 0.5 ml

/// Data Sheet

The Matrix 0,5 ml attachment is suitable for 24 x 0.5 ml reaction vials. The ideal attachment design ensures an optimal temperature transfer and thus guarantees valid results. The thermoshakers Matrix Orbital Delta plus and Matrix Delta plus automatically recognize, which attachment is connected to the device via an interface in the locking mechanism. The attachment can be changed quickly and easily with an ergonomic lever. As part of IKA's Integrated Safety Concept, the insulating plastic coating of the attachment ensures user safety, as high temperatures are not transferred to the frame.

Speed max.: 2000 rpm



designed for scientists

Temperature max.: 100 °C

Suitable for:

Matrix Orbital Delta plus

Matrix Delta plus



designed for scientists

Technical Data

Dimensions (W x H x D) [mm]	141 x 45 x 150
Weight [kg]	0.45