



designed for scientists



## RV 10.4 Dry Ice Condenser

/// Data Sheet

Dry ice condenser for distilling low-boiling solvents. The solvent to be distilled can be continuously fed in through the PTFE inlet tube. Cooling by dry ice, no cooling water required. Max. condensation thanks to low temperatures.



designed for scientists

#### Technical Data

|                                    |                   |
|------------------------------------|-------------------|
| Type of cooling                    | Dry ice condenser |
| Cooling surface [cm <sup>2</sup> ] | 620               |
| Weight [kg]                        | 1.9               |