



designed for scientists



R 1352 Centrifugal stirrer

/// Data Sheet

Two-bladed, blades open with increasing speed. For stirring in round vessels with narrow necks. Effect is similar to that of a 4-bladed propeller stirrer. Medium to high speeds required.

Stirrer Ø: 60 / 15 mm

Shaft Ø: 8 mm

Shaft length: 350 mm



designed for scientists

Technical Data

Shaft diameter [mm]	8
Shaft length [mm]	350
Speed max. [rpm]	2000
Material in contact with medium	AISI 316L, AISI 304
Weight [kg]	0.1362