



designed for scientists



R 4012 SY Propeller stirrer

/// Data Sheet

PTFE-coated propeller stirrer.

Shaft length: 480 mm

Shaft diameter: 10 mm

Dimensions stirring tool (W/H): 50 mm / 10 mm

Suitable for the EasySyn reactor system with reactor vessels SY 1000, SY 1000 D, SY 2000, SY 2000 D (Volumes: 1000 ml; 2000 ml).



designed for scientists

Technical Data

Nominal length [mm]	480
Stirrer design	Pitched blade
Shaft diameter [mm]	10
Drive shaft diameter [mm]	8
Height (agitator unit) [mm]	10
Width (agitator unit) [mm]	50
Depth (agitator unit) [mm]	50
Flow direction	axial
Flow direction for clockwise rotation	upwards
Material in contact with medium	PTFE