




STICKMAX

/// Data Sheet

Universal adhesive mat for shakers.

- Ideal for frequently changing vessel types and sizes
- Self-adhesive
- Vessels can be easily removed thanks to side tilting movement
- Suitable for disinfection

Dimensions (W x D): 200 x 200 mm

Permissible ambient temperature: 5- 80 °C



designed for scientists

Max. speed: 300 rpm



designed for scientists

Technical Data

Dimensions (W x H x D) [mm]	9 x 200 x 200
Weight [kg]	0.136