



SY NS.8 Swivelling screw fitting

/// Data Sheet

Adapter for holding the Temperature sensors PT 100.300 and PT 100.500 and for holding hoses or tubes with an outer diameter of 8 mm.

Fits in NS 29/32 ground joints.

The picked-up object can be fixed in the desired position with an angle of inclination of up to 7°.

Suitable for vacuum applications.

Material: PTFE, PPS.

Temperature resistance: -20 - 200 °C.



designed for scientists

Technical Data

Height (room dimension) [mm]	78
Outer diameter [mm]	40
Vacuum application	yes
Standard ground joint	NS 29/32
Drive shaft diameter [mm]	12
Material in contact with medium	PTFE
Stirring element fastening	chuck